

Enabling High-throughput 3D Cell-based Assays with Commercially-available Sources of Human iPSC-derived Cell Types


Coby Carlson, Rebecca Fiene, Michelle Curtis, David Majewski, Ravi Vaidyanathan, Scott Schachtele, Kirk Twaroski, and Simon Hilcove

FUJIFILM Cellular Dynamics, Inc., Madison, WI USA

FUJIFILM
Value from Innovation

Overview

- Increased assay complexity can be achieved with 3D culture systems.
- Human iPSC technology has helped to enhance physiological relevance by building a bridge between animal testing and human diseases.
- Like “organoids”, 3D spheroid products can be created by mixing individual cell types at defined # and ratios and allowing them to self-assemble in co-culture.
- This approach enables a more flexible method (modular incorporation of defined cell types, including disease-specific lines) while also allowing more control over variability (individual components are well-defined and highly reproducible), both of which are critical to the success of incorporating such technologies into cell-based assay workflows.
- Here we present examples of enabling 3D cell-based assays with human iPSC-derived cell types in 384-well format:
- “Cardiospheres” from cardiomyocytes, cardiac fibroblasts, and endothelial cells show improved myocardial maturity via positive inotropic response (inc. Ca^{2+} waveform amplitude) to compounds like isoproterenol and dobutamine.
- A diverse range of “neurospheres” can be created using healthy or diseased neurons and astrocytes to model neurodegenerative disease, with the option of incorporating microglia to study neuroinflammation.
- Heps and Macs together in 3D yield a more complex liver co-culture system
- Newly released iCell Hepatic Stellate Cells (HSC) are truly quiescent out of thaw and play a critical role in hepatic homeostasis.
- Importantly, all these systems described above are isogenic, meaning the cells are derived from the same iPSC donor background.
- Commercially available iPSC-derived cell types from FCDI are manufactured at scale, quality controlled, and cryopreserved, so that they are ready-to-use with confidence at any point in time.
- Implementing these cells into modular 3D assay workflows is novel approach to complement current organoid research and enhance the biological complexity required for drug discovery, toxicity screening, & disease modeling.




Trending Topics:


- 3D cell culture
- Organoids
- Spheroids
- Advanced models
- Organ-on-a-chip

Materials, Methods, and Instruments


► Induced Pluripotent Stem Cell Technology




► Individual Human iPSC-derived Cell Types = iCell Products




► Ultra Low Attachment (ULA) 3D Cell Culture Spheroid Plates




► Platform Technologies to Monitor, Image, Characterize 3D Spheroids




iCell Cardiospheres: Human, isogenic, tri-culture 3D microtissues for advanced cardiac assays



iCell Neurospheres: Create your own 3D human iPSC-derived “brain-in-a-dish”



iCell Liverspheres: Hepatocytes, Hepatic Stellate Cells, and Macrophages together in 3D



Summary and Future Directions

There is a great deal of excitement and promise around 3D cell culture. Research with iPSC-based organoids and spheroids is focused on creating a more complex and biologically relevant model system to bridge the translational gap, better our understanding of human disease, and facilitate the discovery of new drugs. The advantage of the modular approach presented here with commercially available, cryopreserved, ready-to-use cells is that you can vary cell type numbers and ratios to have exquisite control over the spheroid size and composition. Importantly, this method further enables disease modeling studies and the incorporation of iPSC-derived microglia into neurospheres, for example, is something that is not easily accomplished with organoids. FUJIFILM CDI will continue to develop new applications with iCell products and establish advanced co-culture systems that are amenable to standard workflows.

