

# Development of a Seizure Assessment Model Using Human iPSC-Derived Neural Cells.

○Eriko Watanabe<sup>1)</sup>, Takumi Komikawa<sup>2)</sup>, Yasushi Shiraishi<sup>2)</sup>, Hayato Miyoshi<sup>1)</sup>

1) Bio Science &amp; Engineering Laboratories, FUJIFILM Corporation 2) Process Engineering &amp; Technology Center, FUJIFILM Corporation

## Background

- To design drug compounds without the risk of seizurogenicity is one of the significant issues for drug discovery. Although the risk is usually assessed by animal studies, there are practical difficulties such as the difficulty in assessing seizures through behavioral observations or occurrence of spontaneous seizure.
- In vitro assessment of seizurogenicity using rat primary neural cells or human induced pluripotent stem cell (hiPSC)-derived neural cells have been reported recently. However, verification of the assay is not enough because of the lack of knowledge about species differences, lot-to-lot differences in cells or the effect of variable experimental settings.

## Purpose of research

We established a seizures model using a co-culture system of hiPSC-derived neurons and astrocytes.

This study tested whether the model could detect effects of various compounds with different mechanisms under single or cumulative exposure and whether seizure risks could be identified at clinically toxic concentration ranges.

## Result

### Summary of Data

|      |                                | Single Dose assay |                      |                      |                  | Multiple Dose by Cumulative Addition |                      |                      |                  |               |                      |                      |                  |
|------|--------------------------------|-------------------|----------------------|----------------------|------------------|--------------------------------------|----------------------|----------------------|------------------|---------------|----------------------|----------------------|------------------|
|      |                                | Study A           |                      |                      |                  | Study B                              |                      |                      |                  | Study C       |                      |                      |                  |
|      |                                | Conc. (μg/mL)     | Spikes (times/10min) | Bursts (times/10min) | NB (times/10min) | Conc. (μg/mL)                        | Spikes (times/10min) | Bursts (times/10min) | NB (times/10min) | Conc. (μg/mL) | Spikes (times/10min) | Bursts (times/10min) | NB (times/10min) |
| PTX  | GABA antagonist                | 0                 | 31953                | 265                  | 16               | 0                                    | 19872                | 167                  | 9                | 0             | 73123                | 658                  | 46               |
|      |                                | 10                | 35316                | 393                  | 23               | 0.10                                 | 22832                | 187                  | 10               | 0.316         | 83356                | 678                  | 49               |
|      |                                |                   |                      |                      |                  | 1.1                                  | 22587                | 184                  | 11               | 1             | 88738                | 711                  | 52               |
| PTZ  | GABA antagonist                | 0                 | 19675                | 229                  | 15               | 0                                    | 15215                | 156                  | 9                | 0             | 72848                | 603                  | 40               |
|      |                                | 200               | 19616                | 275                  | 20               | 100                                  | 14746                | 157                  | 10               | 95            | 67046                | 581                  | 38               |
|      |                                |                   |                      |                      |                  | 380                                  | 12461                | 152                  | 10               | 300           | 74769                | 638                  | 42               |
| AXP  | D2R antagonist                 | 0                 | 19675                | 229                  | 15               | 0                                    | 15215                | 156                  | 9                | 0             | 72848                | 603                  | 40               |
|      |                                | 3.0               | 29559                | 947                  | 14               | 0.30                                 | 14746                | 157                  | 10               | 0.316         | 67046                | 581                  | 38               |
|      |                                |                   |                      |                      |                  | 1.3                                  | 17039                | 305                  | 24               | 1             | 101995               | 858                  | 69               |
| DPH  | D2R antagonist                 | 0                 | 19675                | 229                  | 15               | 0                                    | 15215                | 156                  | 9                | 0             | 72848                | 603                  | 40               |
|      |                                | 3.0               | 29559                | 947                  | 14               | 0.30                                 | 14746                | 157                  | 10               | 0.316         | 67046                | 581                  | 38               |
|      |                                |                   |                      |                      |                  | 1.3                                  | 17039                | 305                  | 24               | 1             | 101995               | 858                  | 69               |
| CPZ  | D2R antagonist                 | 0                 | 19675                | 229                  | 15               | 0                                    | 15215                | 156                  | 9                | 0             | 72848                | 603                  | 40               |
|      |                                | 3.0               | 40035                | 960                  | 27               | 0.30                                 | 14746                | 157                  | 10               | 0.316         | 67046                | 581                  | 38               |
|      |                                |                   |                      |                      |                  | 1.3                                  | 17039                | 305                  | 24               | 1             | 101995               | 858                  | 69               |
| LPD  | K <sup>+</sup> channel blocker | 0                 | 19675                | 229                  | 15               | 0                                    | 15215                | 156                  | 9                | 0             | 72848                | 603                  | 40               |
|      |                                | 10                | 15466                | 411                  | 60               | 0                                    | 15215                | 156                  | 9                | 0             | 72848                | 603                  | 40               |
|      |                                |                   |                      |                      |                  | 100                                  | 14746                | 157                  | 10               | 95            | 67046                | 581                  | 38               |
| 4-AP | K <sup>+</sup> channel blocker | 0                 | 19675                | 229                  | 15               | 0                                    | 15215                | 156                  | 9                | 0             | 72848                | 603                  | 40               |
|      |                                | 30                | 17931                | 141                  | 9                | 0                                    | 15215                | 156                  | 9                | 0             | 72848                | 603                  | 40               |
|      |                                |                   |                      |                      |                  | 100                                  | 14746                | 157                  | 10               | 95            | 67046                | 581                  | 38               |
| PILO | M2R agonist                    | 0                 | 19675                | 229                  | 15               | 0                                    | 15215                | 156                  | 9                | 0             | 72848                | 603                  | 40               |
|      |                                | 10                | 19242                | 125                  | 6                | 0                                    | 15215                | 156                  | 9                | 0             | 72848                | 603                  | 40               |
|      |                                |                   |                      |                      |                  | 100                                  | 14746                | 157                  | 10               | 95            | 67046                | 581                  | 38               |
| TPH  | A1R inhibitor                  | 0                 | 19675                | 229                  | 15               | 0                                    | 15215                | 156                  | 9                | 0             | 72848                | 603                  | 40               |
|      |                                | 10                | 19242                | 125                  | 6                | 0                                    | 15215                | 156                  | 9                | 0             | 72848                | 603                  | 40               |
|      |                                |                   |                      |                      |                  | 100                                  | 14746                | 157                  | 10               | 95            | 67046                | 581                  | 38               |
| STRY | GLR antagonist                 | 0                 | 19675                | 229                  | 15               | 0                                    | 15215                | 156                  | 9                | 0             | 72848                | 603                  | 40               |
|      |                                | 10                | 19242                | 125                  | 6                | 0                                    | 15215                | 156                  | 9                | 0             | 72848                | 603                  | 40               |
|      |                                |                   |                      |                      |                  | 100                                  | 14746                | 157                  | 10               | 95            | 67046                | 581                  | 38               |
| AAP  | Antipyretic agent              | 0                 | 19675                | 229                  | 15               | 0                                    | 15215                | 156                  | 9                | 0             | 72848                | 603                  | 40               |
|      |                                | 10                | 19242                | 125                  | 6                | 0                                    | 15215                | 156                  | 9                | 0             | 72848                | 603                  | 40               |
|      |                                |                   |                      |                      |                  | 100                                  | 14746                | 157                  | 10               | 95            | 67046                | 581                  | 38               |

: Significant positive differences between compound and vehicle (DMSO) (Dunnett's test, p&lt;0.05)

: Significant negative differences between compound and vehicle (DMSO) (Dunnett's test, p&lt;0.05)

- ✓ The risk of seizurogenicity was correctly assessed for 8 compounds.
- ✓ The NB was the most sensitive parameter in multiple addition assay.
- ✓ Compounds without clinical data reacted at concentration ranges similar to those reported previously.
- ✓ AAPs that have not been reported to cause seizures clinically were negative.
- ✓ When NB occurs at a high frequency, it is not correctly detected.
- ✓ GABA antagonists (PTX, PTZ) and 4-AP yielded inconsistent results.

## Conclusion

In our developed a seizure model, we detected the effects of test compounds with different mechanisms under single or cumulative exposure. We also confirmed that seizure risks could be identified at concentration ranges reported to cause clinical toxicity. Moving forward, we aim to establish reproducibility.

## Materials and Methods


### Measuring equipment

- ✓ MEA : Maestro Pro (Axion BioSystems)
- ✓ MEA plate: Cytoview (Axion BioSystems)

### Reagents

- ✓ Cell : iCell® GlutaNeurons, iCell® Astrocytes (FUJIFILM Cellular Dynamics, Inc.)
- ✓ Coating agent : 0.1% polyethylenimine
- ✓ Media : BrainPhys Neural Medium with N2, iCell Neural Supplement B, iCell Nervous System Supplement, and Penicillin-Spreptomycin

### Assay Protocol




### Test Compounds

| Compounds               | Mechanism of Action            | Seizure risk | Clinical side effects                     |
|-------------------------|--------------------------------|--------------|---|
| Picrotoxin (PTX)        | GABA antagonist                | Positive     | No data                                   |
| Pentylentetrazole (PTZ) |                                |              | No data                                   |
| Amoxapine (AXP)         |                                |              | No data                                   |
| Diphenhydramine (DPH)   | D2R antagonist                 |              | Toxic range 1.00~4.00 µg/mL <sup>1)</sup> |
| Chlorpromazine (CPZ)    |                                |              | Toxic range 750 µg/mL~ <sup>2)</sup>      |
| Linopiridine (LPD)      | K <sup>+</sup> channel blocker |              | No data                                   |
| 4-Aminopyridine (4-AP)  |                                |              | No data                                   |
| Pilocarpine (PILO)      | M2R agonist                    |              | Lethal dose 60 mg <sup>3)</sup>           |
| Theophylline (TPH)      | A1R inhibitor                  |              | Toxic range 25 µg/mL~ <sup>4)</sup>       |
| Strychnine (STRY)       | GLR antagonist                 |              | No data                                   |
| Acetaminophen (AAP)     | Antipyretic agent              | Negative     | No data                                   |


1. Journal of Japanese Society for Emergency Medicine, 23, 611-615 (2005)  
2. CONFORMIN SUGAR-COATED TABLETS Interview Form Version 14 (2004)  
3. <https://accidentinfo.nih.gov/japanese/mvdsq/13.7.html>  
4. Theophylline SR Capsules 100mg (SANDOZ) Interview Form Version 2 (2004)

### %Change of NetworkBurst

\* : Significant differences between compound and vehicle (DMSO) (Dunnett's test, p&lt;0.05)




### Correlation of Clinical information



- ✓ The Compounds with high brain penetration (DPH, CPZ, and TPH) increase in NB was observed at concentration ranges reported to toxicity range.

### q-PCR



- ✓ mRNA was extracted from the cells four weeks after seeding, and q-PCR was performed.
- ✓ The expression of genes associated with seizurogenicity was confirmed.