(90 × 210mm Size)

**FUJIFILM**

**Wako**

291-77401 ( 1,000 tests)  
297-77403 (10,000 tests)

## Fluorospark™ Kinase/ADP Multi-Assay Kit

### [Introduction]

Phosphorylation and dephosphorylation of proteins in signal transduction play a critical role in physiological regulation as immune response, oncogenesis, differentiation, apoptosis and cell proliferation.

Kinases, which is the enzyme catalyze the phosphorylation of proteins, lipids and other molecules, is closely related to the variety of the diseases including a cancer and to be a main drug discovery target today.

“Fluorospark™ Kinase/ADP Multi-Assay Kit” is a simple and highly sensitive ADP quantitation kit through by enzymatic coupling reaction. The generated ADP through kinase reaction has been quantified as red fluorescent “Resorufin”. By the quantification of Resorufin, this kit can easily quantify the kinase activity with high sensitivity. Furthermore the accuracy of the measured value is higher than previous bioluminescent methods. This kit can utilize for both the endpoint assay and the real-time assay.


### [Features]

- High sensitive assay for kinase activity by quantification of the ADP amount as red fluorescent “Resorufin”.
- Fine linearity until 30 μmol/L ADP.
- Measurement of kinase activity after 30 minutes reaction [Endpoint assay].
- Measurement of kinase activity continuously [Real-time assay].
- Measurement of the wide range of kinases.
- High reliability (Excellent Z'-factor\* at low substrate conversion rate).


※ Z'-factor is one of the static values to evaluate the performance of the HTS assay in terms of sensitivity and accuracy. The Z'-factor > 0.5 indicated the robust assay.

### [Assay Principle]

#### Kinase Reaction

#### Enzymatic Coupling Reaction (ADP Detection Reaction)

- 1/20 -

### [Kit Contents]

No.		1,000 tests*	10,000 tests*
①	Substrate Solution For the enzymatic coupling reaction. [Glucose and NADP]	9mL	90mL
②	Resazurin Solution Resazurin is the precursor of red fluorescent resorufin. Dissolved in DMSO, please make sure the reagents are dissolving well before use.	100 μL	1mL
③	Enzyme Solution For the enzymatic coupling reaction. [ADP Hexokinase (from <i>Thermococcus litoralis</i> ), Glucose-6-phosphate Dehydrogenase (recombinant), Diaphorase I (from <i>Bacillus Stearothermophilis</i> )]	500 μL	5mL
④	Reductant Blocker For inactivation of the reductants (DTT, TCEP etc.). (For the endpoint assay use.) Dissolved in DMSO, please make sure the reagents are dissolving well before use.	400 μL	4mL
⑤	Stop Solution For stopping of the enzymatic coupling reaction. (For the endpoint assay use.)	10mL	100mL
⑥	10mmol/L ATP Solution For the kinase reaction and preparing a standard curve. Please use the provided ATP solution, because the commercial ATP may be contaminated with ADP.	100 μL	1mL
⑦	10mmol/L ADP Solution For preparing a standard curve.	100 μL	1mL

\*For 384 well plate (full volume type).

The test number is depending on the type of microwell plate.

	384 well plate Low volume	384 well plate Full volume	96 well plate
1,000 tests (Code No. 291-77401)	2,000 tests	1,000 tests	200 tests
10,000 tests (Code No. 297-77403)	20,000 tests	10,000 tests	2,000 tests

### [Materials required]

#### < Reagents >

- Ultrapure water  
Avoiding the contamination of ATPase or others.
  - High purity kinase
  - Substrate of kinase
  - The appropriate buffer for kinase reaction  
Phosphate buffer react with stop solution and this causes the fluorescence signal instability. In case of using the phosphate buffer in your kinase reactions, please avoid to use the stop solution.
- The acceptable concentration of the compounds are indicated on the following table.

- 2/20 -

Compounds	Acceptable concentration
NaCl [Sodium Chloride]	400mmol/L
CaCl <sub>2</sub> [Calcium Chloride]	20mmol/L
MgCl <sub>2</sub> [Magnesium Chloride]	100mmol/L
MnCl <sub>2</sub> [Manganese Chloride]	20mmol/L
DMSO [Dimethyl Sulfoxide]	10%
Triton X-100	0.40%
Tween20	0.40%
Na <sub>3</sub> VO <sub>4</sub> [Sodium Orthovanadate(V)]	1mmol/L
NaF [Sodium Fluoride]	20mmol/L
EDTA [Ethylenediaminetetraacetic Acid]	15mmol/L
DTT [Dithiothreitol]	
Endpoint assay	6mmol/L
Real-time assay	1mmol/L

※ Under the condition of 40mmol/L Tris-HCl (pH7.5), 20mmol/L MgCl<sub>2</sub>, 0.1mg/mL BSA.

#### < Instruments >


- Microwell plate**  
The black colored, low-binding microplates is suitable for the fluorescent measurement to avoid the noise signal and the non-specific absorption of proteins and peptides.
- Fluorescence microplate reader**  
For detection of the fluorescence of Resorufin.  
The following microplate readers are available ; Safire, Infinite F200, Infinite F500 [TECAN], EnVision [Perkin Elmer], PHERAstar [BMG LABTECH].  
Filters : Excitation : 540nm, Emission: 590nm (both bandwidth are 10nm)
- Microwell plate mixer**
- Plate seals**

#### 【 Assay methods 】

This kit is used for both the endpoint and the real-time assay.  
[CAUTION] Staurosporine (CAS No. 62996-74-1, serine-threonine kinase inhibitor) inhibits the enzymes including this kit. Please use other kinase inhibitor as the positive control before screening.

#### I . Endpoint Assay

- I-1. Preparation of the mixture of ADP and ATP (for the standard curve)
- I-2. Kinase reaction
- I-3. Preparation of 2 × Detection Reagent
- I-4. Quantitation of ADP amount

This assay method is described for using the low volume 384 well plate. Please adjust the volume of the reagents for other type of microwell plate.

	384 well plate Low volume	384 well plate Full volume	96 well plate
Kinase reaction solution	5 μL	10 μL	50 μL
2 × Detection Reagent	5 μL	10 μL	50 μL
Stop Solution <sup>※</sup>	5 μL	10 μL	50 μL

※ Stop solution : Specific inhibitor for the enzymatic coupling reaction.  
Usually, all of ADP in the enzymatic coupling reaction solution are expended for the resorufin production for 10 minutes after addition of 2 × Detection Reagent. The fluorescence signal of resorufin continues at least 7 hours. If kinase activity has been continued during the enzymatic coupling activity, the fluorescence would be increasing. Addition of stop solution stabilizes the fluorescence.

#### I-1. Preparation of the mixture of ADP and ATP (for the standard curve)

To estimate the ATP-to-ADP conversion rate in your kinase reaction, you can make the standard curve with provided ATP and ADP solutions. In the following steps, you make the ATP/ADP mixture for the standard curve. (Please regard the ATP concentration in your kinase reaction as the 0% conversion rate)

[Example for preparation] ATP concentration ; 10 μmol/L.

- (1) Dilute provided 10mmol/L ADP[⑦] and ATP[⑥] solution.  
Kinase reaction buffer : 10mmol/L ADP (and ATP) = 1,000 : 1
- (2) Mix 10 μmol/L ADP and ATP solution respectively at the rate below.

Conversion rate (%)	0	1	2	3	4	5	10	20	30	40	60	100
10μmol/L ADP(μL)	0	1	2	3	4	5	10	20	30	40	60	100
10μmol/L ATP(μL)	100	99	98	97	96	95	90	80	70	60	40	0
total(μL)	100	100	100	100	100	100	100	100	100	100	100	100

- (3) Add 5 μL of each ADP and ATP mixture solution (the conversion rate is 0~100%) to the low volume 384 well plate (n=2).

#### I-2. Kinase reaction

- (1) Perform kinase reaction with the arbitrary kinase, substrate and reaction buffer and ATP provided.
- (2) Adjust the final volume of kinase reaction mixture to 5 μL in the case of the low volume 384 well plate.

#### I-3. Preparation of 2 × Detection Reagent

- I-3-1. Calculation of the dilution rate of Resazurin Solution  
Please use the calculating formula below in order to suppress the signal of background to achieve high signal/background ratio (S/B ratio).

- ATP concentration at kinase reaction  $< 30 \mu\text{mol/L}$  :  
[Dilution rate of Resazurin Solution] =  $30 / [\text{ATP concentration at kinase reaction} (\mu\text{mol/L})]$
- ATP concentration at kinase reaction  $> 30 \mu\text{mol/L}$  :  
Please use Resazurin Solution without dilution.

Example) Perform kinase reaction with  $10 \mu\text{mol/L}$  ATP.  
[Dilution rate of Resazurin Solution] =  $30 / 10 = 3$   
The dilution rate of Resazurin Solution is threefold.  
It doesn't need to dilute in the case ATP concentration is over  $30 \mu\text{mol/L}$ .

### I-3-2. Preparation of $2 \times$ Detection Reagent

- Thaw the following reagents at room temperature. Substrate Solution[①], Resazurin Solution[②], Enzyme Solution[③] and Reductant Blocker[④]. Resazurin Solution[②] and Reductant Blocker[④] are dissolved in DMSO. Therefore, please make sure these reagents are dissolving well before use.
- Dilute Resazurin Solution[②] by ultrapure water. The dilution rate is calculated at I-3-1.
- Mix these solution ①-④ as the following table.

$2 \times$ Detection Reagent	Ratio (%)
Substrate Solution [①]	90
Diluted Resazurin Solution [②]	1
Enzyme Solution [③]	5
Reductant Blocker [④]	4
Total	100

- Mix well by vortex mixer.

### I-4. Quantification of ADP.

Quantification of ADP in the kinase reaction(I-2.) and the preparation of the standard curve(I-1.) using  $2 \times$  Detection Reagent on step of I-3-2.

- Add  $5 \mu\text{L}$  of  $2 \times$  Detection Reagent into the wells where the kinase reaction was performed and ADP/ATP mixture solution on the low volume 384 well plate.
- Mix well with the microplate mixer.
- Put the plate-seal on the plate and set it in the dark place for 30 minutes at room temperature.  
\*Stop Solution : Add  $5 \mu\text{L}$  of Stop Solution into each wells if you want to stop the enzymatic coupling reaction completely after this reaction.
- Measure the resorfin fluorescence signal with fluorescence microplate reader (Ex 540nm/Em 590nm).


## II. Real-time Assay

The kinase kinetic can be measured with the following steps (real-time assay).

### II-1. Preparation of $2 \times$ RT Detection Reagent

### II-2. Preparation of the kinase substrate solution and the kinase solution with/without kinase

### II-3. Real-time quantification of ADP

(This is the case of assay of using low volume 384 well plate.)

### II-1. Preparation of $2 \times$ RT Detection Reagent

When you perform the real-time assay, please do not dilute the resazurin solution regardless of the ATP concentration in the kinase reaction.

Reductant blocker may inhibit the kinase reactions. Therefore, please do not add the reductant blocker in the  $2 \times$  RT detection reagent.

- Thaw the following reagents at room temperature. Substrate Solution[①], Resazurin Solution[②] and Enzyme Solution[③]. Resazurin is dissolved in DMSO. Therefore, please make sure the reagents is dissolving well before use.
- Mix these solution ①~③ as the following table.

$2 \times$ RT Detection Reagent	Ratio (%)
Substrate Solution [①]	90
Resazurin Solution [②]	1
Enzyme Solution [③]	5
Ultrapure Water	4
Total	100

- Mix well by vortex.

### II-2. Preparation of the kinase and substrate solution with/without kinase

- Preparation of the kinase substrate solution  
Add the kinase substrate and  $10\text{mmol/L}$  ATP solution (attached with this kit) into the kinase reaction buffer.
- Preparation of kinase solution (Assay)  
Dilute the kinase stock solution by the kinase reaction solution.

\* $2 \times$  RT detection reagent do not contain the reductant blocker. Therefore, the reductant such as DTT in the kinase assay solution causes the continuous increase of the background signal during the real-time assay. At least  $1 \text{mmol/L}$  DTT is available for the real-time assay by subtracting the background signal from the assay signal. However, lower reductant concentration should be used for the less background signal.

- (3) Preparation of the Blank solution (Blank : without kinase)  
Prepare the kinase stock solution buffer (without kinase) and dilute it with the kinase reaction solution at the same ratio to (2).

### II-3. Real-time quantification of ADP

- (1) Add 2 × RT Detection Reagent, kinase substrate solution and kinase solution with/without kinases into the low volume 384 well plate.

	Blank [without kinase]	Assay [with kinase]
2 × RT Detection Reagent	5 μL	5 μL
Kinase substrate Solution	2.5 μL	2.5 μL
Kinase solution without kinase	2.5 μL	–
Kinase solution	–	2.5 μL

※This assay is in the case of using low volume 384 well plate. In the case of use other type of microwell plate, please add equal volume of 2 × RT Detection Reagent into combined solution of kinase substrate and kinase solution.

※Blank well is needed to determine the background signal.

- (2) Measure the fluorescence continuously at Ex 540nm/Em 590nm with fluorescence microplate reader.  
(3) Subtract the signal intensity of the blank well from that of the assay well, and calculate the net signal intensity.

※Because it takes about 10 minutes to convert ADP to resorufin in the enzymatic coupling reaction, it is difficult to measure the kinetics of kinases reaction which ends within 10 minutes.

### III. Optimization of ATP concentration, kinase substrate and kinase.

#### III-1. Optimization of ATP concentration

Please adjust the appropriate ATP concentration to get the sufficient signal intensity. In the case of inhibitors screening, high ATP concentration has possibility to interrupt the detection of ATP-competitive inhibitor. Please adjust the ATP concentration as low as in the range to get the sufficient signal intensity.

#### III-2. Optimization of the kinase substrate concentration

Please prepare the dilution series of kinase substrate and set the kinase substrate concentration to get the sufficient signal intensity.

#### III-3. Optimization of kinase concentration

Please prepare the dilution series of kinase and measure the concentration-dependent fluorescent intensity. Please adjust the appropriate kinase concentration which keep the linearity between kinase concentration and fluorescent intensity.

#### [Example for Use]

##### Standard curve of ADP

- Add 5 μL of several concentration ADP solution and 5 μL of 2 × Detection Reagent into the low volume 384 well plate.
- Keep at room temperature for 30 minutes.
- The fluorescent was determined with Safire [TECAN]. The result indicated that the standard curve keep linearity up to 30 μmol/L of ADP.


Fig. 1. Standard curve of ADP

#### ADP standard curve at the low substrate conversion rate.

- Prepare the several standard curve with the substrate conversion rate is 0~30% at ADP and ATP concentration is 1, 10 and 100 μmol/L.
- Add 5 μL of several rate of ADP/ATP solution and 5 μL of 2 × Detection Reagent into the low volume 384 well plate and keep at room temperature for 30 minutes.
- The fluorescent was determined with Safire [TECAN]. The result indicated that the quantitative capability is successful in the condition of the low substrate conversion rate.


Fig. 2. Standard curve of ADP at the low substrate conversion rate.

#### Calculation of Z'-factor

Z'-factor is one of the statical values to evaluate the performance of the HTS assay in terms of sensitivity and accuracy. The Z'-factor at 0% and 5% conversion rate of 1, 10 and 100 μmol/L ADP/ATP concentration was calculated (n=16).

- Add 5 μL of ADP/ATP solution and 5 μL of 2 × Detection Reagent into the low volume 384 well plate and keep at room temperature for 30 minutes.
- The fluorescent was measured with Safire [TECAN]. The result indicated that excellent Z'-factor is acquirable in any ATP concentration.


Fig. 3. Calculation of Z'-factor

### Inhibition curve of PKA by H-89

Preparation of the inhibition curve of cAMP-dependent protein kinase (PKA) by PKA specific inhibitor, H-89.

- 1) Add 5  $\mu\text{L}$  of kinase reaction solution (including PKA 0.02U/ $\mu\text{L}$ , ATP 5  $\mu\text{mol/L}$ , Kempeptide 125  $\mu\text{g/mL}$  and the several concentration of H-89) on the low volume 384 well plate and keep at room temperature for 30 minutes.
- 2) Add 5  $\mu\text{L}$  of 2  $\times$  Detection Reagent and keep at room temperature for 30 minutes.
- 3) The fluorescent was measured with Safire [TECAN]. The result indicated that the determined of  $\text{IC}_{50}$  value was consistent with the reference's value ( $\text{IC}_{50} = 40\text{nmol/L}^{1)}$ .


Fig. 4. PKA inhibition curve by H-89

### Calculation $K_m(\text{ATP})$ of PKA by the real-time assay.

- 1) Add 50  $\mu\text{L}$  of kinase reaction solution (including PKA 0.02U/ $\mu\text{L}$ , the several concentration of ATP and Kempeptide 100  $\mu\text{g/mL}$ ) and 50  $\mu\text{L}$  of 2  $\times$  RT Detection Reagent into the 96 well plate.
  - 2) Measure the fluorescence continuously at room temperature.
- Determine  $K_m(\text{ATP})$  of PKA by calculating the fluorescence increasing per unit time. The result indicated that  $K_m(\text{ATP})$  value of PKA was consistent with the reference's value ( $K_m = 3 \sim 15 \mu\text{mol/L}^{2)}$ .


Fig. 5. Measurement of  $K_m(\text{ATP})$  of PKA

### [References]

- 1) Hidaka, H., Watanabe, M. and Kobayashi, K., "Properties and use of H series compounds as protein kinase inhibitors.", *Methods Enzymol.*, **201**, 328-39(1991).
- 2) Flockhart, D.A. and Corbin, J.D., "Regulatory mechanisms in the control of protein kinases.", *CRC Crit. Rev. Biochem.*, **12**, 133-86(1982).
- 3) Kumagai, K., Kojima, H., Okabe, T. and Nagano, T., "Development of a highly sensitive, high-throughput assay for glycosyltransferase using enzyme-coupled fluorescence detection.", *Anal. Biochem.*, **447**, 146-155(2014).

### [Storage Condition]

Store at  $-20^\circ\text{C}$ . 10mmol/L ATP solution is at risk of decomposing into ADP by repeated freeze-thaw cycle which could be cause to increase its background. At the first thaw, we recommend dispensing into single-use amount to avoid freeze-thaw cycle.

It was checked that the kit components except for 10mmol/L ATP solution have little effect on several times of the freeze/thaw cycle.

## FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://www.wako-chem.co.jp>

### FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

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- (1) Hidaka, H., Watanabe, M. and Kobayashi, K., "Properties and use of H series compounds as protein kinase inhibitors.", *Methods Enzymol.*, **201**, 328-39(1991).
- (2) Flockhart, D.A. and Corbin, J.D., "Regulatory mechanisms in the control of protein kinases.", *CRC Crit. Rev. Biochem.*, **12**, 133-86(1982).
- (3) Kumagai, K., Kojima, H., Okabe, T. and Nagano, T., "Development of a highly sensitive, high-throughput assay for glycosyltransferase using enzyme-coupled fluorescence detection.", *Anal. Biochem.*, **447**, 146-155(2014).

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