


## Myco-Visible Mycoplasma Rapid Test Kit

For Simple, Sensitive, and Rapid Detection of Mycoplasma Contamination in Cell Cultures



**Pack Size:** 24 tests  
**Storage:** -20° C  
**Cat. No.:** 09305090S (5 Tests)  
                  093050901 (24 Tests)  
**Content Version:** Oct 2024

# Table of Contents

1. Introduction to Myco-Visible Mycoplasma Rapid Test Kit .....	3
2. Kit Components and User Supplied Materials .....	4
3. Storage and Kit Stability .....	5
4. Important Consideration Before Use .....	5
5. Protocol.....	5
6. Flow Chart.....	8
7. Data .....	9
8. Troubleshooting .....	11
9. Product Use Limitation & Warranty .....	12

## 1. Introduction to Myco-Visible Mycoplasma Rapid Test Kit

Mycoplasmas are commonly found in research laboratories and are unaffected by many antibiotics used to control bacterial contamination in cell cultures, particularly those targeting cell wall formation, such as penicillin. Their lack of a cell wall also allows some cells to pass through 0.2 µm filters. Several studies have investigated the prevalence of mycoplasma contamination in cell cultures, with estimates ranging as high as 80%. In the USA, contamination rates are lower, around 10-15%. A global average contamination rate between 25% and 50% seems plausible. This significant issue can lead to serious economic losses for the industry. Most cases of cell culture contamination are caused by about six to eight species with human, porcine, or bovine hosts, with *Mycoplasma orale* and *Mycoplasma hyorhinis* being the most common.

Myco-Visible Mycoplasma Rapid Test Kit has been demonstrated to be a high sensitivity, specific and rapid assay for the detection of mycoplasma contamination in cell cultures. Myco-Visible Mycoplasma Rapid Test Kit achieves this by targeting and amplifying unique regions of the mycoplasma genome at loci that are well conserved in all *Mollicutes*, including *Acholeplasma laidlawii*, *Mycoplasma arginini*, *Mycoplasma fermentans*, *Mycoplasma hominis*, *Mycoplasma hyorhinis*, *Mycoplasma orale*, *Mycoplasma salivarium* and *Mycoplasma pneumoniae*. The detection procedure can be completed within 1 hour and the detection limit is as low as 10 Colony-forming Units/mL (CFU/mL) or 10 fg of Mycoplasma genomic DNA per reaction. Specificity test showed no cross-reactivity with other bacterial, fungal and mammalian DNA. The test assay is simple to perform with no special technical expertise or lab equipment. Positive results can be visualized on an immunochromatographic test device in less than 5 minutes after buffer loading.

## 2. Kit Components and User Supplied Materials

### 2.1 Myco-Visible Mycoplasma Rapid Test Kit Components

Components	24 Tests (Cat. No. 093050901)		5 Tests (Cat. No. 09305090S)	
	Package	Cat. No.	Package	Cat. No.
Reaction Beads	3 x 8 ea	093050902	5 ea	09305090S2
Sample Buffer	1 vial (2.5 mL)	093050903	1 vial (0.6 mL)	09305090S3
Primer Mix	1 vial (110 µL)	093050904	1 vial (22 µL)	09305090S4
Nuclease-free Water	1 vial (400 µL)	093050905	1 vial (80 µL)	09305090S5
Positive Control	1 vial (50 µL)	093050906	1 vial (10 µL)	09305090S6
Running Buffer	1 bottle (5 mL)	093050907	1 bottle (1 mL)	09305090S7
Detection Device	24 x 1 ea	093050908	5 x 1 ea	093050908
Quick-start Protocol	1 ea			
Instruction Manual	Available <a href="http://www.mpbio.com">www.mpbio.com</a>			
MSDS & CoA	Available <a href="http://www.mpbio.com">www.mpbio.com</a>			

### 2.2 User Supplied Materials

- Microcentrifuge capable of working at speed of 20,000 x g; if not available, use the maximum speed (>14,000 x g).
- Heat block
- Thermocycler (optional, recommended to be used for 70 °C incubation)
- Sterile, Nuclease-free 0.2 mL and 1.5 mL microcentrifuge tubes
- Pipettes with corresponding filter tips

### 3. Storage and Kit Stability

All Myco-Visible Mycoplasma Rapid Test Kit components are guaranteed until the expiration date printed on the kit label. Upon receipt, store the unopened kit at -20 °C until first use. Once opened, it is recommended to keep the Running Buffer and Detection Device at room temperature for the ease of use. The stability of the two components will not be affected if they are stored at a lower temperature (e.g. 4 °C, -20 °C), but it is required to equilibrate them to room temperature for every use. Always store the other components: Reaction Beads, Primer mix, Nuclease-free Water and Positive Control at -20 °C when not in use.

### 4. Important Consideration Before Use

- ☐ Thaw the reagents at 4 °C or room temperature. Make sure the Detection Device and Running Buffer are at room temperature before use.
- ☐ It is highly recommended to prepare aliquots for Primer Mix (e.g. 5 x 22 µL) and Positive Control (e.g. 5 x 10 µL) upon first use to prevent frequent freeze-thaw and contamination.
- ☐ Clean the work area with 70 % alcohol or 10 % household bleach before the assay setup.
- ☐ Preparation of reaction mix should be performed in a dedicated clean area that is separate from LAMP incubation and sample loading of Detection Device. Do not bring incubated reaction tubes back into the reaction preparation area.
- ☐ Always use filter pipette tips that are sterile and nuclease free.
- ☐ Samples, amplified products and used cassettes are biohazards that must be sterilized before discarding.

## 5. Protocol

### 5.1 Sample Preparation

1. Transfer 1 mL of cell culture supernatant to a 1.5 mL microcentrifuge tube.  
*Note: Test cell culture should be at more than 80% confluency. The supernatant can be stored at 4 °C for up to 3 days or long term at -20 °C before use.*
2. Centrifuge at maximum speed (20,000 x g) for 2 mins and carefully decant the supernatant.  
*Note: The pellet may be invisible. The presence of the pellet will not affect the result.*
3. Add 100 µL of **Sample Buffer** to the tube, mix well by pipetting up and down for 3 times.
4. Heat the sample at 95 °C for 5 min using a heat block.
5. Centrifuge at 14,000 x g for 2 mins and transfer the supernatant to a fresh 1.5 mL centrifuge tube.
6. Use the supernatant immediately as the sample in Section 5.2 or store at -20°C until required.

### 5.2 Reaction Mix preparation

1. Take an appropriate number of **Reaction Beads** provided in individual tubes. The number of **Reaction Beads** includes your samples, a positive control (PC) and a no template control (NTC).
2. Determine the volume of **Nuclease-free water** (X µL) to add to each **Reaction Beads** tube (see table below). Add X µL of **Nuclease-free water** and gently pipette up and down to dissolve the **Reaction Beads**. Add 4 µL of **Primer Mix** and up to 4 µL of sample (or 1 µL control\*) to each tube.

*Note: This step should be performed in a clean and separate location to prevent carry-over contamination.*

Reaction Mix	Volume
Primer Mix	4 µL
Sample / Control*	up to 4 µL / 1 µL
Nuclease-free water	X (to a total of 20 µL)

\* For PC and NTC, use 1 µL of Positive Control and Nuclease-free water respectively (supplied in the kit).

### 5.3 Mycoplasma detection


1. Incubate the reaction mix tube at **70°C** for **30 min** using a thermocycler or heat block.  
*Note: Do not incubate for more than 40 min as it will significantly affect the results.*
2. Allow the reaction tube to equilibrate to room temperature, then briefly spin down the mix.
3. Pipette 5 µL of the reaction mix and load to the sample well of **Detection Device**. Add 3 drops of **Running Buffer** to the sample well and wait for 2 to 5 min for the results.

### 5.4 Result Interpretation


- One band (C-band only): Negative for Mycoplasma.
- Two bands (C-band and T-band): Positive for Mycoplasma.

*Note: A faint T-band is also considered positive for mycoplasma. No control band indicates the test is invalid.*

*Note: The control band will appear within 1-2 min. If positive, the test band will appear within 5 min, test is considered negative if no T-band appeared after 5 min. Any results read after 20 min should be considered invalid since the bands may not be stable after prolonged air exposure.*




## 6. Flow Chart






## 7. Data

Myco-Visible Mycoplasma Rapid Test kit is designed to detect cell culture contaminating Mycoplasma and Acheloplasma species. The detection limit is as low as 10 Colony-forming Units/mL. To identify the analytical sensitivity, inactivated mycoplasma preparations titrated to 100 and 10 CFU/mL from 7 *Mycoplasma* spp in the matrix of a cell suspension of  $1 \times 10^6$  CHO cells/mL were tested in the study. The 7 *Mycoplasma* spp tested here contain key cell culture contaminating and human infecting mycoplasma that are listed in Pharmacopoeia (EP2.6.7/USP <63>/JP G3). The detection limit of the kit (10 CFU/mL) meets the sensitivity guidance of nucleic acid amplification technology (NAT)-based mycoplasma detection method provided in EP/USP/JP.



**Figure 1:** Analytical sensitivity test of Myco-Visible Mycoplasma Rapid Test Kit. All seven strains of mycoplasma were tested at concentrations of 100 and 10 CFU/mL. All samples yielded positive results at both concentrations. The negative control is a no template reaction mix (nuclease-free water) incubated along with other samples.

Myco-Visible Mycoplasma Rapid Test Kit specifically detects mycoplasma without cross-reactivity with other bacterial, fungal or mammalian DNA. In this specificity test, purified gDNA at  $1 \times 10^5$  copies/reaction from various off-target species were assessed. All species give a negative result for mycoplasma detection and there was no interference in the test.



**Figure 2: Specificity test of Myco-Visible Mycoplasma Rapid Test Kit.** In the upper row, purified gDNA at  $1 \times 10^5$  copies/reaction from 10 bacteria species that are genetically related to mycoplasma were tested. In the lower row, gDNA isolated from human, mouse, yeast (15,16) as well as uncontaminated cell culture (13,14) were tested. All samples showed a negative result, indicating no cross-reactivity.

## 8. Troubleshooting

Problem	Possible Cause	Recommendation
<i>Two bands in no template control (NTC)</i>	Reaction time too long	<ul style="list-style-type: none"> <li>Do not exceed incubation time (70 °C for 40 min) as this significantly impacts the reaction.</li> </ul>
	Carryover contamination of the reagent	<ul style="list-style-type: none"> <li>Discard current aliquot of reagent and use new one.</li> <li>Reaction set up should be performed in dedicated environments free of contamination and separate from the incubation and result interpretation site.</li> <li>Use of filter tips will reduce the contamination.</li> </ul>
<i>Only Control band is present in Positive Control</i>	Degradation of reagent	<ul style="list-style-type: none"> <li>Avoid frequent freeze-thaw cycles of the Primer Mix and Positive Control.</li> </ul>
	Weak <i>Mycoplasma</i> infection of the sample	<ul style="list-style-type: none"> <li>Increase the volume of template (up to 5 µL).</li> <li>Grow the cell culture for additional 48 hours and repeat the test.</li> <li>Use first-time LAMP product as a template to repeat the test.</li> </ul>
<i>No bands</i>	Detection Device spoiled	<ul style="list-style-type: none"> <li>Store the device under proper condition.</li> <li>Equilibrate the device to room temperature before use.</li> <li>Use within 1 hour after opening</li> </ul>
	Running Buffer contamination	<ul style="list-style-type: none"> <li>Ensure the Running Buffer is stored properly and not contaminated.</li> </ul>
<i>Weak Test band</i>	Mild contamination	<ul style="list-style-type: none"> <li>If a light pink test band is observed for test sample and no band observed for NTC. The results should be considered positive.</li> <li>If both the test sample and NTC show weak test bands, the result should be considered invalid due to contamination.</li> </ul>

## 9. Product Use Limitation & Warranty

The products presented in this instruction manual are for research or manufacturing use only. They are not to be used as drugs or medical devices in order to diagnose, cure, mitigate, treat or prevent diseases in humans or animals, either as part of an accepted course of therapy or in experimental clinical investigation. These products are not to be used as food, food additives or general household items. Purchase of MP Biomedicals products does not grant rights to reproduce, modify, or repackage the products or any derivative thereof to third parties. MP Biomedicals makes no warranty of any kind, expressed or implied, including merchantability or fitness for any particular purpose, except that the products sold will meet our specifications at the time of delivery.

Buyer's exclusive remedy and the sole liability of MP Biomedicals hereunder shall be limited to, at our discretion, no replacement or compensation, product credits, refund of the purchase price of, or the replacement of materials that do not meet our specification. By acceptance of the product, Buyer indemnifies and holds MP Biomedicals harmless against, and assumes all liability for, the consequence of its use or misuse by the Buyer, its employees, or others, including, but not limited to, the cost of handling. Said refund or replacement is conditioned on Buyer notifying within thirty (30) days of receipt of product. Failure of Buyer to give said notice within thirty (30) days shall constitute a waiver by the Buyer of all claims hereunder with respect to said material(s).

**Australia**

Tel: +61 2.8824.2100  
Tel: +61 1800.249.998  
Email: [custserv.au@mpbio.com](mailto:custserv.au@mpbio.com)

**Austria & Germany**

Tel: 0800.426.67.337  
Tel: 00800.7777.9999  
Email: [custserv.de@mpbio.com](mailto:custserv.de@mpbio.com)

**Belgium**

Tel: 00800.7777.9999  
Email: [custserv.be@mpbio.com](mailto:custserv.be@mpbio.com)

**Canada**

Tel: +1 800.854.0530  
Email: [custserv.ca@mpbio.com](mailto:custserv.ca@mpbio.com)

**China**

Tel: +86 400.150.0680  
Email: [custserv.cn@mpbio.com](mailto:custserv.cn@mpbio.com)

**Europe**

Tel: +33 3.88.67.54.25  
Tel: +33 00800.7777.9999  
Email: [custserv.eur@mpbio.com](mailto:custserv.eur@mpbio.com)

**France**

Tel: +33 3.88.67.54.25  
Email: [custserv.fr@mpbio.com](mailto:custserv.fr@mpbio.com)

**India**

Tel: +91 22.27636921/22/25  
Email: [custserv.in@mpbio.com](mailto:custserv.in@mpbio.com)

**Italy**

Tel: 00800.7777.9999  
Email: [custserv.it@mpbio.com](mailto:custserv.it@mpbio.com)

**Japan**

Tel: +81 3.6667.0730  
Email: [custserv.jp@mpbio.com](mailto:custserv.jp@mpbio.com)

**Latin America**

Tel: +1 800.854.0530  
Tel: +1 440.337.1200  
Email: [custserv.la@mpbio.com](mailto:custserv.la@mpbio.com)

**New Zealand**

Tel: +64 9.912.2460  
Email: [custserv.nz@mpbio.com](mailto:custserv.nz@mpbio.com)

**North America**

Tel: +1 800.854.0530  
Tel: +1 440.337.1200  
Email: [custserv.na@mpbio.com](mailto:custserv.na@mpbio.com)

**Poland**

Tel: 00800.7777.9999  
Email: [custserv.po@mpbio.com](mailto:custserv.po@mpbio.com)

**Russia**

Tel: +7 495 604.13.44  
Email: [custserv.rs@mpbio.com](mailto:custserv.rs@mpbio.com)

**Serbia**

Tel: +381 11.242.1972  
Email: [custserv.se@mpbio.com](mailto:custserv.se@mpbio.com)

**Singapore/ APAC**

Tel: +65 6775.0008  
Tel: +65 6394.7675  
Email: [custserv.ap@mpbio.com](mailto:custserv.ap@mpbio.com)

**South Korea**

Tel: +82 2.425.5991  
Email: [custserv.kr@mpbio.com](mailto:custserv.kr@mpbio.com)

**Switzerland**

Tel: 00800.7777.9999  
Email: [custserv.ch@mpbio.com](mailto:custserv.ch@mpbio.com)

**The Netherlands**

Tel: 00800.7777.9999  
Email: [custserv.nl@mpbio.com](mailto:custserv.nl@mpbio.com)

**United Kingdom**

Tel: 0800.282.474  
Email: [custserv.uk@mpbio.com](mailto:custserv.uk@mpbio.com)

[www.mpbio.com](http://www.mpbio.com)

