


ヒト近位尿細管上皮細胞三次元培養モデルを用いた薬物誘発性腎毒性評価

Etsushi Takahashi¹, Kaoru Morimura¹, Ayano Araki¹, Daichi Higuchi², Hiroshi Arakawa², Ikumi Tamai², Yoichi Jimbo¹¹ R & D Department, Precision Engineering Center, Industrial Division, Nikkiso Co. Ltd² Faculty of Pharmaceutical Sciences, Institute of Medical Pharmaceutical and Health Sciences, Kanazawa University

Introduction


Kidney-derived cell lines¹⁾ or Primary human kidney cells²⁾ do not maintain the gene expression related to kidney function such as drug transporters³⁾. Conventional kidney cells have not yet been used in drug discovery. Therefore, many of the drug-induced kidney injury (DIKI) has been evaluated by animal studies. However, *in vitro* evaluation of DIKI using human cells is desired from the viewpoint of low predictivity to the clinical trial, species difference and animal welfare.

We investigated the usefulness of evaluating DIKI using three-dimensional cultured human proximal tubular epithelial cells (3D-RPTEC), whose expression levels of major drug transporters are comparable to those of human kidney cortex. To detect for DIKI more sensitively, we investigated the usefulness of High Content Analysis (HCA) using a confocal image cytometer.



Methods

2D-culture : RPTEC (LONZA, CC-2553, Passage 3) was thawed and cultured in REGM (LONZA).
3D-culture : RPTECs were cultured as spheroids in ultra-low attachment 96-well culture plate (PrimeSurface, Sumitomo Bakelite). The medium was changed once every 2-3 days. This spheroid of RPTEC (3D-RPTEC®) will be available from Nikkiso from July 2023.



Intracellular ATP content was measured using the CellTiter-Glo® 3D Cell Viability Assay (Promega) by addition to the wells of a 96-well plate after drug exposure of 3-28 days. For High Content Analysis (HCA), the images of cells were taken using a confocal image cytometer CQ1 (Yokogawa).

Conclusion

● 新規なヒト腎細胞3D-RPTEC®を用いた*in vitro*腎毒性モデル

- 長期培養できる腎細胞 (薬物トランスポーターやメガリンを発現)
- トランスポーターを介した腎毒性を検出

● ATP測定によりヒト腎毒性リスクを評価

- ATP測定とHCAを組み合わせて高い確度で腎毒性を検出

問い合わせ先

3D-RPTEC®について興味またはご質問のある方は以下までお問い合わせ下さい。

日機株式会社

<https://www.nikkiso.co.jp/products/industrial/3drptec/>




メールアドレス

3D-RPTEC@nikkiso.co.jp

Results


Fig.1 Features of 3D-RPTEC

(A) Drug transporter (Microarray)




3D-RPTEC showed improved expression of drug transporters and was maintained in long-term culture (Fig.1-A).

(B) Endocytosis (Proteomics)



Proteomics analysis of plasma membranes showed higher expression of megalin and cubilin, which are responsible for endocytosis in kidney, compared to 2D-culture (Fig.1-B).

(C) Donor-to-donor variability



(D) Immunocytochemistry






Fig.3 Safety margin for DIKI

(A) ATP assay in 3D-RPTEC



(B) ATP assay in HK-2



The safety margin of DIKI was calculated from the concentration that inhibits ATP by 20% (IC₂₀ value) and the effective blood concentration (C_{max}).

Sensitivity: 66.7% (day7), 76.2% (day28)


Specificity: 100% (day7), 100% (day28)

Compound	C _{max} total (μM)	7 days exposure IC ₂₀ (μM)	IC ₂₀ / C _{max}
Mitomycin C	1.77	0.24	0.14
Carboplatin	67.3	14.5	0.22
Oxaliplatin	4.37	1.35	0.31
Cisplatin	7.47	2.93	0.39
Cephalothin	54.7	127	2.32
Rifampicin	12.2	35.4	2.91
Pentamidine	0.9	3.7	4.11
Cephaloridine	192	812	4.22
Zoledronic acid	0.97	10.5	10.8
Vancocycin	38	1140	30
Tenofovir	0.93	44.2	47.68
Gentamycin	15.7	878	55.9
Omeprazole	3.18	212	66.5
Tetracycline	3.67	256	69.9
Doxorubicin	0.04	4.8	>100
Cyclosporin A	0.21	>100	>100
Amphotericin B	1.52	>100	>100
Ifosfamide	199	>300	>100
Adefovir	0.07	77.6	>100
Tacrolimus	0.03	>10	>100
Allopurinol	73.5	>3000	>100
Phenacetin	38.4	>3000	>100
Ciprofloxacin	12.1	>1000	>100
Lisinopril	0.09	>10	>100
Isoproterenol HCL	0.95	>100	>100
Mannitol	43.1	>10000	>100
Metformin	12.4	>1000	>100
Thiamine	6.8	>1000	>100
Dexamethasone	0.68	>100	>100
Dopamine	0.08	>10	>100
Atropine	0.09	>10	>100
Pindolol	0.2	>30	>100

* : Schulz, M. et al., Crit Care.(2012)16: R136

Fig.4 High Content Analysis

(A) Target molecule in DIKI



(B) Image analysis






Fig.5 DIKI evaluation by combination of ATP and HCA

(A) Evaluation method in HCA



(B) Combination of HCA and ATP



* : $p < 0.05$ (chi-square test vs. ROC curve AUC=0.500)

COI Disclosure Information

I have no COI regarding this presentation.
Presenting author: Etsushi Takahashi
Organization: Nikkiso Co. Ltd

Acknowledgement

We gratefully thank Naoki Ishiguro (Nippon Boehringer Ingelheim Co) and Nagisa Kato (Yokogawa Electric Corporation) for supporting the study. This research was supported by Japan Agency for Medical Research and Development (AMED).

References

(1) Jenkinson et al., Eur J Physiol. (2012) 464:601-611
(2) C. Welland et al., Toxicology in Vitro (2007)21:466-491