

Testing for *BRCA1* and *BRCA2* Gene Mutations

***BRCA* Gene Mutations and Hereditary Cancer Risk**

Breast cancer is the most common type of cancer in women after skin cancer; in the general population, 1 in 8 women (about 13%) will develop breast cancer in their lifetime. While many factors may increase the risk of cancer, such as smoking or exposure to chemicals or radiation, about 5-10% of cancer cases are associated with inherited genetic mutations ([National Cancer Institute, 2019](#)). Genes called ***BRCA1*** and ***BRCA2*** (short for BReast CAncer gene) were the first genes to be associated with a strong family history of breast and ovarian cancer.




Figure 1. Chromosomal locations of *BRCA1* and *BRCA2* genes.
(Tessssa13, CC BY-SA 4.0, via Wikimedia Commons.)

Everyone has the *BRCA1* and *BRCA2* genes (Figure 1), which code for proteins that normally serve important functions in repairing damaged DNA. But women with certain mutations in the *BRCA1* gene have up to a 72% lifetime risk of breast cancer and 44% risk of ovarian cancer (compared to a 1% risk of ovarian cancer for the general population), while women with *BRCA2* mutations have up to a 69% lifetime risk of breast cancer and 17% risk of ovarian cancer ([National Cancer Institute, 2020](#)). *BRCA2* mutations also increase the risk of breast and prostate cancer in men. Individuals with *BRCA* gene mutations tend to be diagnosed with more aggressive cancer at a younger age, often in their 30s and 40s.

Genetic Testing and Personalized Medicine

Because people with *BRCA* gene mutations have such significantly increased risks of breast and ovarian cancer, testing for *BRCA* gene mutations can be an important step in personalizing cancer screening and risk management based on genetic risk factors. Individuals with a strong family history of breast and/or ovarian cancer may be referred by their doctor to a genetic

counselor, who will take a detailed family medical history and use mathematical models to estimate that individual's personal cancer risk. Based on the results of this cancer risk assessment, the genetic counselor may recommend testing for mutations in the *BRCA* genes and possibly other additional genes associated with increased risk of breast and ovarian cancer. Women who test positive for *BRCA* gene mutations typically undergo more frequent screening for breast cancer, usually involving both a mammogram and an MRI every year. Some women with *BRCA* gene mutations and a strong family history of cancer decide to undergo preventive mastectomies and oophorectomies (removal of ovaries) to reduce their future risk of developing cancer. Angelina Jolie notably raised public awareness about risk-reducing surgery when she discussed her personal decision to have a preventive mastectomy given her *BRCA1* mutation status ([Jolie, 2013](#)).

The Promise of Next Generation Sequencing

To test for *BRCA* gene mutations, a patient sample, typically a blood sample, is sent to a lab that specializes in genetic testing. The goal of genetic testing is to examine the DNA sequence of the gene to determine if any harmful mutations are present. Traditionally, this has involved amplifying many copies of the gene of interest (*BRCA1* or *BRCA2* in this case) using the polymerase chain reaction (PCR), followed by Sanger sequencing using fluorescently-labeled nucleotides to determine the DNA sequence. While Sanger sequencing is still considered the gold standard to confirm DNA mutations, it can be time-consuming to sequence large genes. The length of a typical Sanger sequencing read is approximately 1 kb, while the *BRCA1* gene is 5.6 kb in length and *BRCA2* is 10.3 kb ([Nicolussi et al., 2019](#)).

Recent advances in DNA sequencing technology have allowed scientists to utilize **next generation sequencing (NGS)** as an alternative approach to traditional Sanger sequencing methods for detecting genetic mutations ([Nicolussi et al., 2019](#), Figure 2). While Sanger sequencing reads a single DNA sequence at a time, an NGS approach first breaks the DNA into many smaller fragments (200-500 bp), to which specific identifier tags are added at the ends to create a library of DNA fragments from the gene of interest. Each DNA library uses unique tags that tie the sequences back to specific genes from a specific patient. Millions of DNA fragments from multiple sample libraries (e.g. multiple patients being tested for both *BRCA1* and *BRCA2* mutations) can then be amplified and sequenced at a time using a high-throughput NGS platform. The DNA sequences of overlapping fragments are then aligned using software to assemble the DNA sequence of the entire gene. Because NGS can sequence many DNA samples simultaneously, using NGS technology for genetic testing can increase the number of different genes and patient samples that can be tested while decreasing turnaround time and sequencing cost per sample, enabling doctors and patients to receive *BRCA* genetic testing results more quickly and cost effectively.




Figure 2. Workflow for analyzing *BRCA1* and *BRCA2* mutations using NGS.

DNA is extracted from a blood sample, then libraries of tagged DNA fragments are prepared for NGS. After amplification and sequencing, DNA analysis includes alignment of overlapping DNA fragments and identification of genetic variants (mutations). Sanger sequencing is used both to troubleshoot regions that did not amplify well and to confirm variants identified by NGS. A technique called MLPA is used to identify large genomic rearrangements that can be missed by NGS.

([Nicolussi et al., 2019](#)).

BioChain Launching New *BRCA*-Characterized Breast Cancer Tissue Lines

Stay tuned for BioChain's upcoming blog on its new line of *BRCA1*- and *BRCA2*-characterized breast cancer tissue samples! These tissue samples have been characterized for *BRCA1* and *BRCA2* genetic mutations using NGS, and can be used in genetic testing validation, diagnostic assay development, cellular localization studies, and many other applications.

BioChain offers many other [breast cancer-related products](#), including tumor tissue arrays with normal tissue controls, triple negative breast cancer tissue samples, and CancerSeq tumor tissue lines prescreened by NGS for single nucleotide polymorphisms, insertions, deletions and copy number variation. BioChain also offers [custom procurement and characterization services](#) to meet specific research needs.

References:

1. National Cancer Institute. (2019, March 15). The Genetics of Cancer. <https://www.cancer.gov/about-cancer/causes-prevention/genetics>
2. National Cancer Institute. (2020, November 19). BRCA Gene Mutations: Cancer Risk and Genetic Testing. <https://www.cancer.gov/about-cancer/causes-prevention/genetics/brca-fact-sheet>
3. Jolie, A. (2013, May 14). My Medical Choice. *The New York Times*, A25. <https://www.nytimes.com/2013/05/14/opinion/my-medical-choice.html>
4. Nicolussi, A., Belardinilli, F., Mahdavian, Y., Colicchia, V., D'Inzeo, S., Petroni, M., Zani, M., Ferraro, S., Valentini, V., Ottini, L., Giannini, G., Capalbo, C., & Coppa, A. (2019). Next-generation sequencing of *BRCA1* and *BRCA2* genes for rapid detection of germline mutations in hereditary breast/ovarian cancer. *PeerJ*, 7, e6661. [10.7717/peerj.6661](https://doi.org/10.7717/peerj.6661).



Dr. Lianna Wong received her Ph.D. in Molecular and Cell Biology from the University of California, Berkeley. She has 20 years of experience in biology education at Santa Clara University and UC Berkeley.

Metadata:

HTML Description:

Scientists can use next generation sequencing (NGS) to more efficiently detect genetic mutations in genes like *BRCA1* and *BRCA2*.

HTML Keywords: *BRCA1*, *BRCA2*, breast cancer, personalized medicine, genetic testing, next generation sequencing, NGS, BioChain