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1. Corvari V, Narhi LO, Spitznagel TM, et al. Subvisible (2-100 μm) particle analysis during biotherapeutic drug product development: Part 2, experience with the application of subvisible particle analysis. *Biologicals*. 2015; 43(6): 457-473. doi:10.1016/j.biologicals.2015.07.011
2. Roesch A, Zölls S, Stadler D, et al. Particles in Biopharmaceutical Formulations, Part 2: An Update on Analytical Techniques and Applications for Therapeutic Proteins, Viruses, Vaccines and Cells. *J Pharm Sci*. 2021;000. doi:10.1016/j.xphs.2021.12.011
3. Simon CG, Borgos SE, Calzolari L, et al. Orthogonal and complementary measurements of properties of drug products containing nanomaterials. *Journal of Controlled Release*. 2023;354:120-127. doi:10.1016/j.jconrel.2022.12.049
4. Carpenter JF, Randolph TW, Jiskoot W, et al. Overlooking Subvisible Particles in Therapeutic Protein Products: Gaps That May Compromise Product Quality. *J Pharm Sci*. 2009;98:1201-1205. doi:10.1002/jps
5. Ripple DC, Hu Z. Correcting the Relative Bias of Light Obscuration and Flow Imaging Particle Counters. *Pharm Res*. 2016; 33(3): 653-672. doi:10.1007/s11095-015-1817-9
6. United States Pharmacopeia. USP <787> Subvisible Particulate Matter in Therapeutic Protein Injections.
7. United States Pharmacopeia. USP<788> Particulate Matter in injections.
8. Zölls S, Gregoritz M, Tantipolphan R, et al. How Subvisible Particles Become Invisible — Relevance of the Refractive Index for Protein Particle Analysis. *Pharm Biotechnol*. 2013;102(5):1434-1446. doi:10.1002/jps
9. Shibata H, Harazono A, Kiyoshi M, Ishii-Watabe A. Quantitative Evaluation of Insoluble Particulate Matters in Therapeutic Protein Injections Using Light Obscuration and Flow Imaging Methods. *J Pharm Sci*. 2021;000. doi:10.1016/j.xphs.2021.09.047
10. Strehl R, Rombach-Riegraf V, Diez M, et al. Discrimination between silicone oil droplets and protein aggregates in biopharmaceuticals: A novel multiparametric image filter for sub-visible particles in microflow imaging analysis. *Pharm Res*. 2012;29(2):594-602. doi:10.1007/s11095-011-0590-7

11. Calderon CP, Daniels AL, Randolph TW. Deep Convolutional Neural Network Analysis of Flow Imaging Microscopy Data to Classify Subvisible Particles in Protein Formulations. *J Pharm Sci*. 2018;107(4):999-1008. doi:10.1016/j.xphs.2017.12.008
12. Shibata H, Terabe M, Shibano Y, et al. A Collaborative Study on the Classification of Silicone Oil Droplets and Protein Particles Using Flow Imaging Method. *J Pharm Sci*. 2022;111(10):2745-2757. doi:10.1016/j.xphs.2022.07.006
13. United States Pharmacopeia. USP <1788> Methods for the Determination of Subvisible Particulate Matter.
14. Yokogawa Fluid Imaging Technologies. Measuring Subvisible Particles and Aggregates using FlowCam LO [White Paper]. Published online 2021. <https://info.fluidimaging.com/flowcam-lo-white-paper-download>
15. Yokogawa Fluid Imaging Technologies. Characterizing Protein Aggregates with Orthogonal and Complementary Analytical Techniques [White Paper]. Published online 2022. <https://info.fluidimaging.com/protein-aggregate-orthogonal-methods-white-paper>
16. Narhi LO, Corvari V, Ripple DC, et al. Subvisible (2-100 μm) particle analysis during biotherapeutic drug product development: Part 1, considerations and strategy. *J Pharm Sci*. 2015;104(6):1899-1908. doi:10.1002/jps.24437
17. Telikepalli SN, Carrier MJ, Ripple DC, et al. An Interlaboratory Study to Identify Potential Visible Protein-Like Particle Standards. *AAPS PharmSciTech*. 2023;24(1). doi:10.1208/s12249-022-02457-9

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Copyright © Yokogawa Electric Corporation