(90 × 210mm Size)

FUJIFILM

Wako

Code No. 299-78301 (96 tests)

Human ES/iPS Cell Monitoring Kit

rBC2LCN is a recombinant protein of N-terminal domain of BC2L-C lectin derived from *Burkholderia cenocepacia*. rBC2LCN exhibits significant affinity to a mucin-type O-glycan sugar chain called H-type3 (Fuc α 1-2Gal β 1-3GalNAc) displayed on the podocalyxin binds to the surface of human pluripotent stem cells (hPSCs), human ES cells and human iPS cells. Therefore rBC2LCN is reported as a marker of hPSCs.

The glycoprotein recognized by rBC2LCN is released from hPSCs into cell culture supernatant. This kit can quantitatively determine the released glycoprotein by rBC2LCN-antibody sandwich assay and estimate the number of hPSCs. Because, cell culture supernatants, but not cells itself, are used as sample targeted for assay, it is easily capable of monitoring the increase and decrease of hPSCs in a noninvasive manner during continuous cell culture. This kit is useful for researches, such as regenerative medicine, cell differentiation and iPS cell induction.

This product is for laboratory use only ; use in any such application is the responsibility of the user.

Kit Contents

	Code No.	Contents	Amount
(1)	296-78311	rBC2LCN-coated Plate	96 wells × 1 plate (8 wells × 12 strips)
(2)	293-78321	HRP-conjugated Antibody Solution	400 μ L × 1 tube
(3)	290-78331	Negative Control	1 mL × 1 tube
(4)	297-78341	Positive Control	100 μ L × 1 tube
(5)	294-78351	Dilution Solution	15 mL × 1 bottle
(6)	291-78361	Wash Solution (10 ×)	100 mL × 1 bottle
(7)	298-78371	TMB Solution	10 mL × 1 bottle
(8)	295-78381	Stop Solution	10 mL × 1 bottle
(9)	292-78391	Plate Seal	4 sheets

Storage

Store at 2-10 °C

Package

Code No.	Packaging
299-78301	96 tests

Additional required apparatus

- Centrifuge (possible at 1,700 × g (about 3,000 rpm))
- Vortex mixer
- Plate mixer (desirable)
- Automated 96 well-plate washer (desirable)
- 96 well-plate reader (absorbance 450 nm, reference wavelength 600-650 nm)


Precautions

< About assay >

- Harvest the cell culture supernatants on the next day after the replacement of the culture medium and, at the same time, count the number of the cells in the dish. When the medium in the dish is 5 mL and it is found to be 5×10^6 cells/dish, for example, the number of cells in the cell culture supernatant is regarded to be 1×10^6 cells/mL.
- Generate a standard curve by cell culture supernatants under the undifferentiated culture of hPSCs, and calculate the number of hPSCs in the sample based on the standard curve. However, the calculated number of hPSCs is not always equal to the actual cell number of hPSCs. Please use as one of indicators to monitor the undifferentiated state of cells and the progress of cell differentiation.
- The relation between signal intensity and cell number (cells/mL) is greatly affected by culture conditions such as cell strains and culture medium*. Prepare standard curves for each cell strain and undifferentiated maintenance culture condition, respectively.

※ *The comparison of different culture medium is indicated below.*

After culturing human iPS cells in the Medium A, B, C and D under undifferentiated condition, the culture supernatants were harvested on the day following the date of medium replacement. The cell number of each condition was counted, and the culture supernatant was diluted with individual culture medium to 2,000 cells/mL. The signal intensity of different culture supernatant was measured with this kit.

< About use of kit >

- Wear gloves and laboratory coats, and use laboratory glassware when handling this kit. In particular, carefully handle Stop Solution which is strong acid. If the reagents of the kit touch skin directly, wash with water. When getting them in your eye, don't rub, and please wash with running water and seek treatment immediately.
- Don't use kit past expiration date.
- Don't mix the reagents from other lot or sources.
- Store at 20-25 °C before use.
- The plate seal is sticky. When peeling off the plate seal, please be careful for the solution in the well not to jump up and spill.
- After washed with an automated washer and so on, invert the plate and blot it against clean paper towels to remove residual washing buffer in the wells.
- Store TMB solution in the dark.
- Store kit contents except for the plate seal at 2-10 °C immediately after use.
- Dispose reagents in accordance with applicable regional/national regulatory requirements.

< Preparation of "sample" (cell culture supernatant) >

- Harvest cell culture supernatants cultured for 18-24 hours after medium exchange.
 - ※ Because the culturing time affect the signal intensity, exchange whole culture medium. And harvest cell culture supernatants in constant time.
- Centrifuge the cell culture supernatants at 1,700 × g (about 3,000 rpm) at room temperature for 10 minutes. Transfer the supernatants to new tubes and regard them as measuring "samples".
- When not measuring the "samples" immediately, store at -20 °C.

< Preparation of standard curve >

- Prepare a standard curve for each cell strain and undifferentiated maintenance culture condition.
- Culture hPSCs for 18-24 hours after changing the whole culture medium with the fresh medium and then harvest the cell culture supernatant. Then, detach the hPSCs from the dish and count the cell number of hPSCs. Prepare the "samples" for standard curve from cell culture supernatants according to Preparation of "sample".
- Store the "samples" at -20 °C prior to use (several times of freezing-thawing cycle are possible).
- Dilute "sample" for standard curve with the same new medium with the sample targeted for assay, and make a standard curve. Prepare by three-fold serial dilutions from 30,000 cells/mL to about 41 cells/mL and make sure of the tendency. After that, produce a dilution series by two-fold serial dilutions from the appropriate cell number to make standard curve. Measure only culture medium because culture medium show high background.
- In each assay, use one strip for standard curve, and calculate the cell number of hPSCs in "sample" targeted for assay from the standard curve.

- When it's difficult to use one of strips for making a standard curve because of multiple "samples", decrease the number of wells used for the standard curve (2-4 wells). But, select the number of cells with good linearity of the standard curve.

Procedures

< Preparation of assay >

- Use rBC2LCN-coated plate, Negative Control, Positive Control, Dilution Solution, Wash Solution (10 ×), Stop Solution and Plate Seal after returning to room temperature (except for HRP-conjugated Antibody Solution and TMB Solution).
- Dilute the Wash Solution (10 ×) 10 times with purified water (or deionized water) of room temperature. Wash Solution (1 ×) need at least 40 mL per 1-strip. Prepare the required amount of Wash Solution (1 ×) according to the number of "sample". When using an automated plate washer, prepare an amount of Wash Solution (1 ×) sufficient for washer setting.
- The precipitate occasionally appears in Wash Solution (10 ×). After the precipitate is dissolved by warming to 25-37 °C, please dilute Wash Solution (10 ×).
- Just prior to use, dilute the Positive Control 40 times with the Dilution Solution. Positive Control is used to make sure of success of the reaction. (Add 5 μL of Positive Control to 195 μL of Dilution Solution, mix by vortex mixer and add 50 μL/well to blank wells.)
- Add the Negative Control to a well without diluting.
- Just prior to use, dilute the HRP-conjugated Antibody Solution 20 times with the Dilute Solution. Diluted HRP-conjugated Antibody Solution need 50 μL/well. After transferring the requirement to a tube. Store the rest of HRP-conjugated Antibody Solution at 2-10 °C immediately.
- Before 20-30 minutes of coloring reaction, transfer the required amount of TMB Solution to a sterilized tube, and avoid light exposure, store at room temperature before use. Store the rest of TMB Solution at 2-10 °C immediately.

< Assay >

- (1) After rBC2LCN-coated Plate was returned to the room temperature, take the plate from the bag. Remove unused strips from the plate frame, return them to the bag, and store at 2-10 °C.
- (2) Close the holder of the plate frame and fix the strips. Add 350 μL/well of Wash Solution (1 ×) to each well and aspirate. Repeat the wash process 3 times, then invert the plate and blot it against clean paper towels to remove residual washing buffer in wells.
- (3) Add 50 μL/well of "sample" for standard curve, "sample" targeted for assay, Negative Control and diluted Positive Control to each well. After mixing lightly by a plate mixer, seal the plate with a Plate Seal and incubate for 1 hour at room temperature.
- (4) Remove the Plate Seal, add with 350 μL/well of Wash Solution (1 ×) to each well and aspirate. Repeat the wash process 3 times, and then invert plate and blot it against clean paper towels to remove residual washing buffers in wells.

- (5) Add 50 μ L/well of diluted HRP-conjugated Antibody Solution to each well. After mixing lightly by a plate mixer, seal the plate with a Plate Seal and incubate for 1 hour at room temperature.
- (6) Remove the Plate Seal, add 350 μ L/well of Wash Solution (1 \times) to each well and aspirate. Repeat the wash process 6 times, and then invert the plate and blot it against clean paper towels to remove residual washing buffers in wells.
- (7) Add 50 μ L/well of TMB Solution to each well. After mixing lightly by a plate mixer and incubate for 30 minutes at room temperature (cover it with an aluminum foil in order to protect from light).
- (8) Add 50 μ L/well of Stop Solution to each well and mix lightly to stop the reaction. Read the absorbance of each well at 450 nm (reference wavelength 600-650 nm) by a plate reader. Ensure no bubbles are present in the wells prior to reading plate.
- (9) Calculate the cell number of hPSCs contained in "sample" targeted for assay based on the standard curve.
 ※ When using Positive Control and Negative Control, confirm that the absorbance of Positive Control is more than 0.5 and the absorbance of Negative Control is less than 0.15.

Flow chart

rBC2LCN-coated Plate

↓ Wash plate 3 times

Add 50 μ L of "samples" for standard curve, "sample" targeted for assay, Negative control and diluted Positive Control

↓ Mix and incubate at room temperature for 1 hour

↓ Wash plate 3 times

Add 50 μ L of HRP-conjugated Antibody Solution diluted 20 times

↓ Mix and incubate at room temperature for 1 hour

↓ Wash 6 times

Add 50 μ L of TMB Solution

↓ Mix and incubate at room temperature for 30 minutes (protect from light)

Add 50 μ L of Stop Solution

↓ Mix

Read the absorbance of 450 nm (reference wavelength 600-650 nm)

Reference

Tateno, H., Onuma, Y., Ito, Y., Hiemori, K., Aiki, Y., Shimizu, M., Higuchi, K., Fukuda, M., Warashina, M., Honda, S., Asashima, M. and Hirabayashi, J. : *Sci. Rep.*, **4**, 4069 (2014).

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
 Telephone : + 81-6-6203-3741
 Facsimile : + 81-6-6201-5964
<http://www.wako-chem.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
 Richmond, VA 23237
 U.S.A.
 Telephone : + 1-804-271-7677
 Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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