


Reagent for VOC Analysis

Volatile organic compounds (VOC) are organic chemicals that have a high vapor pressure at ordinary room temperature. VOC are widely used because they are important substances as solvents and fuels. However, when released into the environment, it causes health hazards such as pollution. So VOCs are regulated by law in some countries. We provide a variety of mixture standard solution and columns to analyze VOC.


Mixture Standard Solution of VOC

We provide 4 types of VOC mixture standard solutions. Especially, 25 VOC Mixture Standard Solution is traceable to Japanese measurement standards (Japan Calibration Service System; JCSS ⇒ next page).

■ GC Analysis (25 VOC Mixture Standard Solution)



Column	DB-624 3.0 μ m 0.53mm \times 75m		
Column temp.	30 $^{\circ}$ C(40min) \rightarrow 2 $^{\circ}$ C/min \rightarrow 50 $^{\circ}$ C \rightarrow 5 $^{\circ}$ C/min \rightarrow 60 $^{\circ}$ C(18min) \rightarrow 5 $^{\circ}$ C/min \rightarrow 90 $^{\circ}$ C(10min) \rightarrow 5 $^{\circ}$ C/min \rightarrow 150 $^{\circ}$ C(4min)		
Sample	2 μ L		
Carrier gas	He 4.9mL/min		
Injection	cool on-column	Detector	FID, 250 $^{\circ}$ C
Instrument	Agilent 7890B		



Column	DB-WAX 1.0 μ m 0.53mm \times 60m		
Column temp.	40 $^{\circ}$ C(2min) \rightarrow 20 $^{\circ}$ C/min \rightarrow 80 $^{\circ}$ C(16min) \rightarrow 10 $^{\circ}$ C/min \rightarrow 180 $^{\circ}$ C		
Sample	2 μ L		
Carrier gas	He 2.9mL/min (20min) \rightarrow 4.6mL/min (Flow rate : 1.0mL/min)		
Injection	cool on-column	Detector	FID, 230 $^{\circ}$ C
Instrument	Agilent 7890B		

Peak No.	Constituent	Peak No.	Constituent	Peak No.	Constituent
1	1,1-Dichloroethylene	10	1,2-Dichloroethane	19	<i>o</i> -Xylene
2	Dichloromethane	11	Trichloroethylene	20	Bromoform
3	<i>trans</i> -1,2-Dichloroethylene	12	1,2-Dichloropropane	21	<i>p</i> -Dichlorobenzene
4	<i>t</i> -Butyl Methyl Ether	13	Bromodichloromethane	22	1,4-Dioxane
5	<i>cis</i> -1,2-Dichloroethylene	14	<i>cis</i> -1,3-Dichloropropene	23	<i>p</i> -Xylene
6	Chloroform	15	Toluene	24	<i>m</i> -Xylene
7	1,1,1-Trichloroethane	16	<i>trans</i> -1,3-Dichloropropene	25	1,1,2-Trichloroethane
8	Carbon Tetrachloride	17	Tetrachloroethylene		
9	Benzene	18	Dibromochloromethane		

Product List

Code No.	Product Name	Grade	Volume
229-01651	8 VOC Mixture Standard Solution (each 1mg/mL Methanol Solution)	for Environment Analysis	2mL \times 5A
228-01481	12 VOC Mixture Standard Solution (each 1mg/mL Methanol Solution)	for Environment Analysis	2mL \times 10A
220-02421	14 VOC Mixture Standard Solution (each 1mg/mL Methanol Solution)	for Environment Analysis	2mL \times 5A
225-02351	25VOC Mixture Standard Solution (Methanol Solution)	JCSS	2mL \times 5A

■ Composition List


Constituent	for Environment Analysis			JCSS
	8 VOC 229-01651	14 VOC 220-02421	12 VOC 228-01481	25 VOC 225-02351
Tetradecane	○			
Styrene	○			
Ethylbenzene	○			
Chloroethylene		○		
Carbon Tetrachloride		○	○	○
1,2-Dichloroethane		○	○	○
1,1-Dichloroethylene		○	○	○
<i>trans</i> -1,2-Dichloroethylene		○		○
<i>cis</i> -1,2-Dichloroethylene		○	○	○
1,3-Dichloropropene*		○	○	○
Dichloromethane		○	○	○
Tetrachloroethylene		○	○	○
1,1,1-Trichloroethane		○	○	○
1,1,2-Trichloroethane		○	○	○
Trichloroethylene		○	○	○
Benzene		○	○	○
1,4-Dioxane				○
Chloroform				○
Dibromochloromethane				○
Bromodichloromethane				○
Bromoform				○
Toluene	○			○
<i>t</i> -Butyl Methyl Ether				○
<i>o</i> -Xylene	○			○
<i>m</i> -Xylene	○			○
<i>p</i> -Xylene	○			○
1,2-Dichloropropane				○
<i>p</i> -Dichlorobenzene	○			○

* mixture of isomers

What is JCSS ?

The JCSS reference material is traceable to the specified SI traceable reference material, which is the national measurement standard of Japan, manufactured by the designated calibration institute "Chemicals Evaluation and Research Institute : CERI" under the Measurement Law, and has a JCSS symbol by the JCSS registered business operator. It is provided to the user as a reference material, which attached a certificate with JCSS symbol. We have been certified by IA Japan as an international MRA accreditation operator, and we have issued a calibration certificate with the JCSS certification symbol. The values described in this calibration can be accepted internationally through the MRA of ILAC/APAC.

■ Traceability of JCSS



for Stinking Substances Analysis

Stinking substances cause unpleasant odors and may damage the living environment. In Japan, 22 stinking substances have been designated by the Environment Agency, and all factories within the regulated area designated by the each local government are required to measure odor compounds to considerate of the surrounding environment.

■ Measurement Method of Stinking Substances

No.	Stinking Substances	Regulation Value [ppm]	Analytical Method
1	Ammonia	1-5	Absorption Spectrophotometry
2	Methyl Mercaptan	0.02-0.01	GC
3	Hydrogen Sulfide	0.02-0.2	
4	Dimethyl Sulfide	0.01-0.2	
5	Dimethyl Disulfide	0.009-0.1	GC
6	Trimethylamine	0.005-0.07	
7	Acetaldehyde	0.05-0.5	GC, GC/MS, HPLC
8	Propionaldehyde	0.05-0.5	
9	<i>n</i> -butyraldehyde	0.009-0.08	
10	Isobutyraldehyde	0.02-0.2	
11	<i>n</i> -valeraldehyde	0.009-0.05	
12	Isovaleraldehyde	0.003-0.01	
13	Isobutanol	0.9-20	GC, GC/MS
14	Ethyl Acetate	3-20	
15	Methyl Isobutyl Ketone	1-6	
16	Toluene	10-60	
17	Styrene	0.4-2	
18	Xylene	1-5	
19	Propionic Acid	0.03-0.2	GC
20	<i>n</i> -Butyric Acid	0.001- 0.006	
21	<i>n</i> -Valeric Acid	0.0009-0.004	
22	Isovaleric Acid	0.001-0.01	

Product List


Code No.	Product Name	Grade	Volume
012-08241	Ammonia Standard Solution (0.2μg/μL Water Solution)	Ref. for Stinking Substances Analysis	5mL
012-20081	Ammonia Standard Solution II (100μg/μL Water Solution)	Ref. for Stinking Substances Analysis	5mL
130-06173	Methyl Mercaptan Standard Solution (1μg/μL Benzene Solution)	Ref. for Stinking Substances Analysis	2mL×5A
137-06183	Dimethyl Sulfide Standard Solution (0.1μg/μL Benzene Solution)	Ref. for Stinking Substances Analysis	2mL×5A
040-17253	Dimethyl Disulfide Standard Solution (0.1μg/μL Benzene Solution)	Ref. for Stinking Substances Analysis	2mL×5A
200-06483	Trimethylamine Standard Solution (1μg/μL Ethanol Solution)	Ref. for Stinking Substances Analysis	2mL×5A
012-17391	2 Aldehydes-DNPH Mixed Standard Solution (each 0.1μg/μL Acetonitrile Solution)	Ref. for Stinking Substances Analysis (HPLC)	2mL×5A
018-17491	2 Aldehydes-DNPH Mixture Standard Solution (each 0.1μg/μL Ethyl Acetate Solution)	Ref. for Stinking Substances Analysis (GC)	2mL×5A
012-15451	6 Aldehydes-DNPH Mixture Standard Solution (each 0.1μg/μL Acetonitrile Solution)	Ref. for Stinking Substances Analysis (HPLC)	2mL×5A
012-14851	6 Aldehydes Mixture Standard Solution (each 0.1μg/μL Ethyl Acetate Solution)	Ref. for Stinking Substances Analysis (GC)	2mL×5A
019-27811	6 Aldehydes-DNPH Mixture Standard Solution II (each 0.1mg/mL Acetonitrile Solution)	Ref. for Stinking Substances Analysis (HPLC)	2mL×5A
018-18231	16 Aldehydes-DNPH Mixed Standard Solution (each 10μg/mL Acetonitrile Solution)	Ref. for HPLC	2mL×5A
153-01861	8 Offensive Odor Organic Solvents Mixture Standard Solution each 1μg/μL Pentane Solution) (Isobutanol, Ethyl Acetate, Methyl Isobutyl Ketone, Toluene, Styrene, <i>o</i> -, <i>m</i> -, <i>p</i> -Xylene)	Ref. for Stinking Substances Analysis (GC)	2mL×5A
205-12901	Toluene Standard Solution (1μg/μL Pentane Solution)	Ref. for Stinking Substances Analysis (GC)	2mL×5A
197-10751	Styrene Standard Solution (1μg/μL Pentane Solution)	Ref. for Stinking Substances Analysis (GC)	2mL×5A
246-00661	Xylene Standard Solution (<i>o</i> -, <i>m</i> -, <i>p</i> -each 1μg/μL Pentane Solution)	Ref. for Stinking Substances Analysis (GC)	2mL×5A
168-15141	Propionic Acid Standard Solution (1μg/μL Water Solution)	Ref. for Stinking Substances Analysis	2mL×5A
029-10881	<i>n</i> -Butyric Acid Standard Solution (1μg/μL Water Solution)	Ref. for Stinking Substances Analysis	2mL×5A
221-01111	<i>n</i> -Valeric Acid Standard Solution (1μg/μL Water Solution)	Ref. for Stinking Substances Analysis	2mL×5A
093-03431	Isovaleric Acid Standard Solution (1μg/μL Water Solution)	Ref. for Stinking Substances Analysis	2mL×5A

■ Composition List (Aldehyde Mixture Standard Solution)

Constituent	2 Aldehydes		6 Aldehydes		16 Aldehydes
	012-17391, 018-17491		012-15451	012-14851, 019-27811	018-18231
Formaldehyde	○	○		○	
Acetaldehyde	○	○	○	○	
Propionaldehyde		○	○	○	
Acrolein				○	
Acetone				○	
<i>n</i> -Butyraldehyde			○	○	
Isobutyraldehyde		○	○	○	
Crotonaldehyde				○	
<i>n</i> -Valeraldehyde		○	○	○	
Isovaleraldehyde		○	○	○	
Benzaldehyde				○	
<i>n</i> -Hexanal				○	
<i>o</i> -Tolualdehyde				○	
<i>m</i> -Tolualdehyde				○	
<i>p</i> -Tolualdehyde				○	
2,5-Dimethylbenzaldehyde				○	

■ HPLC Column (Wakopak® Wakosil-DNPH)

Wakopak® Wakosil DNPH specializes in the analysis of aldehydes. This column can simultaneously analyze *iso*-, *n*-butyraldehyde (difficult-to-separate) and other aldehydes. In addition, by using the Wakosil® DNPH Eluent, it is possible to suppress baseline fluctuations and elution of unknown peaks.




- | | |
|-------------------------------------|---------------------------------------|
| 1. Formaldehyde-2,4-DNPH | 10. <i>n</i> -Valeraldehyde-2,4-DNPH |
| 2. Acetaldehyde-2,4-DNPH | 11. Benzaldehyde-2,4-DNPH |
| 3. Propionaldehyde-2,4-DNPH | 12. Hexaldehyde-2,4-DNPH |
| 4. Acrolein-2,4-DNPH | 13. <i>o</i> -Tolualdehyde-2,4-DNPH |
| 5. Acetone-2,4-DNPH | 14. <i>m</i> -Tolualdehyde-2,4-DNPH |
| 6. IsoButyraldehyde-2,4-DNPH | 15. <i>p</i> -Tolualdehyde-2,4-DNPH |
| 7. <i>n</i> -Butyraldehyde-2,4-DNPH | 16. 2,5-Dimethylbenzaldehyde-2,4-DNPH |
| 8. Crotonaldehyde-2,4-DNPH | |
| 9. IsoValeraldehyde-2,4-DNPH | |

Analysis Condition

Column	Wakopak® Wakosil-DNPH (4.6 mmφ x 250 mm)								
Eluent	A) Wakosil® -DNPH Eluent A B) Wakosil® -DNPH Eluent B								
Gradient	<table border="1"> <tr> <th>Time (min)</th> <th>B conc. (%)</th> </tr> <tr> <td>0-16</td> <td>10</td> </tr> <tr> <td>16-35</td> <td>10-90</td> </tr> <tr> <td>35-40</td> <td>90</td> </tr> </table>	Time (min)	B conc. (%)	0-16	10	16-35	10-90	35-40	90
Time (min)	B conc. (%)								
0-16	10								
16-35	10-90								
35-40	90								
Flow Rate	0.8mL/min. at 35 °C								
Inj. Vol.	10 μL								
Detector	UV 360 nm								

■ SPE Column (Presep®-C DNPH)

Presep® DNPH is a SPE column which collects aldehydes and derivatives them with 2,4-dinitrophenylhydrazine (DNPH). This column can reduce noise for lower aldehyde (formaldehyde, acetaldehyde, and acetone).



- | | |
|--|---------------------------------------|
| 1. Formaldehyde-2,4-DNPH | 9. Isovaleraldehyde-2,4-DNPH |
| 2. Acetaldehyde-2,4-DNPH | 10. <i>n</i> -Valeraldehyde-2,4-DNPH |
| 3. Acetone-2,4-DNPH | 11. <i>o</i> -Tolualdehyde-2,4-DNPH |
| 4. Acrolein-2,4-DNPH | 12. <i>m</i> -Tolualdehyde-2,4-DNPH |
| 5. Propionaldehyde-2,4-DNPH | 13. <i>p</i> -Tolualdehyde-2,4-DNPH |
| 6. Crotonaldehyde-2,4-DNPH | 14. Hexaldehyde-2,4-DNPH |
| 7. <i>n</i> -Isobutyraldehyde-2,4-DNPH | 15. 2,5-dimethylbenzaldehyde-2,4-DNPH |
| 8. Benzaldehyde-2,4-DNPH | |

Analysis Condition

Column	Wakopak® Wakosil-DNPH-II (4.6 mmφ x 150 mm)								
Eluent	A) Wakosil® -DNPH-II Eluent A B) Wakosil® -DNPH-II Eluent B								
Gradient	<table border="1"> <tr> <th>Time (min)</th> <th>B conc. (%)</th> </tr> <tr> <td>0-10</td> <td>0</td> </tr> <tr> <td>10-10.01</td> <td>0-100</td> </tr> <tr> <td>10.01-16</td> <td>100</td> </tr> </table>	Time (min)	B conc. (%)	0-10	0	10-10.01	0-100	10.01-16	100
Time (min)	B conc. (%)								
0-10	0								
10-10.01	0-100								
10.01-16	100								
Flow Rate	1.0 mL/min. at 35 °C								
Inj. Vol.	10 μL								
Detector	UV 360 nm								

Product List

Code No.	Product Name	Grade	Volume
238-59411	Wakopak Wakosil-DNPH 4.6*250mm	-	1 column
237-61733	Wakopak Wakosil-DNPH-II 4.6*150mm	-	1 column
290-34251	Presep-C DNPH	Ref for Sample Pretreatment	20 pcs
291-43951	Presep-C DNPH (Short)	Ref for Sample Pretreatment	20 pcs
233-01611	Wakosil® DNPH Eluent A	for HPLC	1 L
230-01621	Wakosil® DNPH Eluent B	for HPLC	1 L
236-02181	Wakosil® DNPH-II Eluent A	for HPLC	1 L
233-02191	Wakosil® DNPH-II Eluent B	for HPLC	1 L

Listed products are intended for laboratory research use only, and not to be used for drug, food or human use. / Please visit FUJIFILM Wako Laboratory Chemicals site: <https://labchem-wako.fujifilm.com/> / This leaflet may contain products that cannot be exported to your country due to regulations. / Bulk quote requests for some products are welcomed. Please contact us.

FUJIFILM Wako Laboratory Chemicals site
<https://labchem-wako.fujifilm.com>



FUJIFILM Wako Pure Chemical Corporation
 1-2, Doshomachi 3-Chome, Chuo-ku, Osaka 540-8605, Japan
 Tel: +81 6 6203 3741 Fax: +81 6 6203 1999
 ffwk-cservice@fujifilm.com

FUJIFILM Wako Chemicals U.S.A. Corporation
 1600 Bellwood Road, Richmond, VA 23237, U.S.A.
 Toll-Free (U.S. only): +1 877 714 1920
 Tel: +1 804 271 7677 Fax: +1 804 271 7791
 wkuslabchem@fujifilm.com

FUJIFILM Wako Chemicals Europe GmbH
 Fuggerstr 12, 41468 Neuss, Germany
 Tel: +49 2131 311 0 Fax: +49 2131 311 100
 labchem_wkeu@fujifilm.com

FUJIFILM Wako Chemicals (Hong Kong) Limited
 Room 1111, 11/F, International Trade Centre, 11-19 Sha Tsui Road,
 Tsuen Wan, N.T., Hong Kong
 Tel: +852-2799-9019 Fax: +852-2799-9808
 wkhk.info@fujifilm.com

FUJIFILM Wako (Guangzhou) Trading Corporation
 Room 3003, 30/F., Dong Shan Plaza 69, Xian Lie Zhong Road,
 Guangzhou, 510095, China
 Tel: +86-20-8732-6381(Guangzhou) Tel: +86-21-6288-4751(Shanghai)
 Tel: +86-10-6413-6388(Beijing)
 wkgz.info@fujifilm.com