FUJIFILM**Wako**

Code No. 296-80901

for Skin Sensitisation Test ADRA Kit

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ADRA (Amino acid Derivative Reactivity Assay) is in chemico test method, which is an alternative (non-animal) method, for predicting skin sensitization potential of chemicals, using *N*-(2-(1-naphthyl)acetyl)-L-cysteine (hereinafter referred to as NAC) and α -*N*-(2-(1-naphthyl)acetyl)-L-lysine (hereinafter referred to as NAL). This reagent kit enables simple implementation of the ADRA.

【Reagent supplied】

This kit contains reagents to implement the ADRA for two 96-well plates.

Components	Amount
NAC (297-80931)	for 10mL × 2 vials
NAL (294-80941)	for 10mL × 2 vials
NAC Buffer (pH 8.0) (premixed) (293-80911)	for 300mL × 2 vials
NAL Buffer (pH 10.2) (premixed) (290-80921)	for 300mL × 2 vials
0.01mol/L EDTA Solution (291-80951)	1mL × 2 vials

【Reagents required other than the kit】

- 1) Trifluoroacetic Acid (TFA) (Wako Special Grade 204-02743)
- 2) Acetonitrile (for HPLC 019-08631)
- 3) Distilled water or higher-grade water*
*Use water with low metal content for the ADRA.
- 4) Acetone** (Guaranteed Reagent 016-00346)
- 5) Dimethyl Sulfoxide (DMSO)** (Guaranteed Reagent 043-07216)
**Not to be used if the test chemical dissolves in water or acetonitrile.
- 6) Phenylacetaldehyde (CAS No. 122-78-1)

【Necessary apparatus】

- 1) Electronic balance –with a displayed precision of 0.1 mg
- 2) Micropipette – Three micropipettes capable of distribution of 2-10 μ L, 10-100 μ L, 100-1000 μ L
- 3) 12 channel pipette – capable of distribution of 50-150 μ L
- 4) HPLC system – equipped with a light-shielding auto-sampler for a 96-well plate capable of liquid feeding at 0.3 mL/min
- 5) UV detector-photodiode array (PDA) detector or absorbance detector (281 nm)
- 6) HPLC column
- 7) pH meter – with precision of ± 0.01 pH, equipped with a buffer solution for calibration
- 8) Incubator – capable of a temperature of 25°C
- 9) 96-well plate
- 10) 500 mL plastic bottle
- 11) Test tube mixer
- 12) Plate seal*
*Use the seal having high sealability and solvent-resistant performance.
- 13) Plate shaker

—1/16—



14) Plate centrifuge

※ All the containers and instrument (except for components of HPLC) used in this test must be disposable products made of Polypropylene or Polyethylene. If contaminated by metal ion, dimerization of the NAC occurs and the assay cannot be performed correctly.

【Principle】

This method assays the presence of skin sensitization by reacting the test chemicals with NAC and NAL for 24 hours at 25°C, then determining the concentration of residual NAC and NAL in the reaction liquid.

NAC and NAL are compounds consisting of a naphthalene ring at a detection area and cysteine or lysine at the nucleophilic area. The relative concentration of NAC and NAL after the 24-hour reaction is determined by UV detection at 281 nm under the conditions of gradient elution in HPLC.

**【Pre-work: Solubility Assessment of Test Chemicals】**

Prior to performing the ADRA, the solubility of the test chemical must be evaluated.

The following 4 solvents are suitable for the solubility evaluation and the priority for the selection of the solvent is as follows.

1. Water*, 2. Acetonitrile, 3. Acetone, 4. DMSO**

*In case of acid anhydride, due to its hydrolysis reaction, water cannot be selected as a solvent.

**DMSO is only used for the 20 mM solution preparation.

Procedure for Solubility Assessment

- 1) Calculate the quantity of the test chemical necessary to prepare 5 mL of 20 mM solution of the test chemical. Use the formula below to determine the amount.

【Preparation of 20 mM solution (5 mL)】

$$\text{Quantity of the test chemical (mg)} = \text{Molecular weight} \times \frac{100}{\text{purity (\%)}} \times \frac{20\text{mM}}{1000} \times 5\text{mL}$$

• If there is no information on the purity of the test chemical, assume that it is 100%.

• If the determined quantity is less than 10 mg, weigh more than 10 mg.

- 2) Weigh the quantity of the test chemical directly in a 10–15 mL centrifuge tube, and record the actual weighed value.

- 3) Calculate the necessary amount of the solvent from the weighed value using the following formula, then prepare the solution.

【Formula to determine the solvent amount for preparation of 20 mM solution】

$$\text{Solvent amount (mL)} = \text{Amount of weighted test chemical (mg)} \times \frac{\text{Purity (\%)}}{\text{Molecular weight}} \times 0.5$$

- 4) If the test chemical is difficult to dissolve, perform ultrasonication for 5 minutes.

- 5) Check visually whether the test chemical is soluble or not.

※ If the test chemical at concentration of 20 mM does not dissolve in all the solvent, perform the solubility check at concentration of

—2/16—

1 mM with water, acetonitrile and acetone. DMSO is only used as solvent for the chemical at concentration of 20mM, as it needs to be diluted with acetonitrile for the reaction test.

【Formula to determine the quantity of the test chemical to prepare 100 mL of 1 mM solution】

$$\text{Quantity of the test chemical (mg)} = \text{Molecular weight} \times \frac{100}{\text{purity (\%)}} \times \frac{1\text{mM}}{1000} \times 100\text{mL}$$

If the test chemical is an aqueous solution, use water as a test solvent, and determine the necessary amount of the test chemical and water using the formula below. If the specific gravity is not available, assume it as 1.

【Formula to determine necessary amount for preparation of 5 mL of 20 mM solution of test chemical】

$$\frac{\text{Molecular weight}}{\text{Concentration (w/v\%)}} \times \frac{10}{\text{Specific gravity}} = \text{Amount of test chemical (\mu L)}$$

$$5\text{ mL} - \frac{\text{Amount of test chemical (\mu L)}}{1000} = \text{Solvent amount (mL)}$$

※ In case of preparing 1 mM solution, dilute 20 mM solution by 20 times.

【Preparation of reagents】

1. Preparation of test chemical solution

Prepare the test chemical solution immediately before the reaction process.

- Prepare the test chemical solution at concentration of 20 mM or 1 mM, using the solvent selected in "Solubility Assessment of Test Chemicals". For preparation of the test chemical solution, refer to the preparation method described in "Procedure for Solubility Assessment".
- In case of 20 mM test chemical solution, when using a solvent other than DMSO, prepare the concentration to be 1 mM by diluting with same solvent by 20 times. When using DMSO as a solvent, dilute with acetonitrile by 20 times to prepare the concentration to be 1 mM.

※ If sediment or cloudiness is observed, stop the experiment at that point and perform the Solubility Assessment of Test Chemical again.

2. Preparation of buffer

- ① Preparation method of 100 mM NAC Buffer (pH 8.0)
 - 1) Measure 300 mL of water.
 - 2) Transfer the content of the NAC Buffer (pH 8.0) premixed into a 500 mL plastic bottle.
 - 3) Take 30 mL from the measured water, and rinse the container of the NAC Buffer (pH 8.0) premixed with it, then pour into the 500 mL plastic bottle.
 - 4) Add 270 mL of the remaining water into the 500 mL plastic bottle, and dissolve the premix.
 - 5) After the premix is dissolved, add 10 μL of 0.01 mol/L EDTA solution.
 - 6) Transfer some of the prepared buffer to another container and check its pH. Ensure that the pH is within the range of 7.9-8.1.

② Preparation method of 100 mM NAL Buffer (pH 10.2)

- 1) Measure 300 mL of water.
- 2) Transfer the content of the NAL Buffer (pH 10.2) premixed into a 500 mL plastic bottle.
- 3) Take 30 mL from the measured water, and rinse the container of the NAL Buffer (pH 10.2) premixed with this, then pour into the 500 mL plastic bottle.
- 4) Add 270 mL of the remaining water into the 500 mL plastic bottle, and dissolve the premix.
- 5) Transfer some of the prepared buffer to another container and check its pH. Ensure that the pH is within 10.1-10.3.
※ If storing the prepared buffer, store it in the refrigerator and use within 3 days after preparation.

3. Preparation of Positive Control Solution

Prepare 1 mM Phenylacetaldehyde (Molecular weight = 120.15) acetonitrile solution as positive control.

• Preparation procedure

- 1) Using the following formula, calculate the quantity of the positive control necessary to prepare 5 mL of 20 mM positive control solution.

$$\text{Molecular weight} \times \frac{10}{\text{purity (\%)}} = \text{Quantity of Phenylacetaldehyde (mg)}$$

- 2) Weigh the quantity of Phenylacetaldehyde directly in a 10-15 mL centrifuge tube.
- 3) Calculate the necessary amount of the solvent based on the actual weighted value, using the following formula.

$$\text{Weighed value of Phenylacetaldehyde (mg)} \times \frac{\text{Purity (\%)}}{\text{Molecule weight}} \times 0.5 = \text{Solvent amount (mL)}$$

- 4) Dissolve the weighted Phenylacetaldehyde in the determined amount of acetonitrile. (concentration: 20 mM)
- 5) Dilute it with acetonitrile by 20 times in an appropriate container* such as 1.4 mL sample track. (concentration: 1 mM)
*96-well deep well plate can also be used instead of sample track.

4. Preparation of NAC solution and NAL solution

Prepare the NAC solution and NAL solution immediately before the reaction process.

- 1) Preparation of NAC solution
Take 10 mL of the NAC Buffer (pH 8.0) and directly add to the container of NAC, then dissolve (concentration: 6.667 μM).
※ Do not perform ultrasonication and intense stirring as it can cause oxidization of NAC.
- 2) Preparation of NAL solution
Take 10 mL of the NAL Buffer (pH 10.2) and directly add to the container of NAL, then dissolve (concentration: 6.667 μM).

5. Preparation of stop solution

Add 1 mL of TFA to 40 mL of water. (concentration: 2.5% (v/v))

【Reaction process】

The following operations are performed in a 96-well plate.

1. Implementation of stability test (0 hour)
A stability test (0 hour) should be implemented to ensure the stability of NAC and NAL.

Refer to the plate layout and dispense the reagents as shown below.

	ST (0)
NAC Solution / NAL Solution	150 μ L
Acetonitrile	50 μ L
Stop Solution (2.5% (v/v) TFA Solution)	50 μ L

After dispensing the reagents, seal the plate with a plate seal, and mix it on a plate shaker. After being centrifuged, implement HPLC analysis immediately.

※ Check the ADRA-specialized webpage for the conditions of HPLC analysis and analysis method for obtained data. The URL is shown on the last page of this manual.

2. Dispensation of reagents

① Refer to the plate layout and dispense the reagents to the corresponding wells as shown below.

(The plate layout which can be used for a 96-well plate is available on the ADRA-specialized webpage.)

	Sample	RC-A	RC-B	RC-C	CC	PC	ST(24)
NAC Solution / NAL Solution	150 μ L	150 μ L	150 μ L	150 μ L	—	150 μ L	150 μ L
NAC Buffer (pH 8.0)/ NAL Buffer (pH 10.2)	—	—	—	—	150 μ L	—	—
Acetonitrile	—	50 μ L	50 μ L	—	—	—	50 μ L
Solvent for dissolving test chemical	—	—	—	50 μ L	—	—	—
Test chemical solution	50 μ L	—	—	—	50 μ L	—	—
Phenylacetaldehyde solution	—	—	—	—	—	50 μ L	—

After the dispensation is complete, seal the plate with a plate seal, and mix it on a plate shaker.

Following spinning down in a centrifuge, incubate it under the light-shielding condition for 24 hours (\pm 1 hour) at 25 $^{\circ}$ C (\pm 1 $^{\circ}$ C).

② Reaction stop

After the 24-hour reaction, add 50 μ L of the stop solution to each well. The stop solution should be added within 30 minutes after the reaction is finished.

Following the adding of the stop solution, seal the plate with a plate seal, and stir it on a plate shaker.

*If precipitation is observed in the reaction liquid after the stirring, centrifuge the plate at a low speed (100-400 g) to gather the deposit at the bottom. If it is possible to collect more than 100 μ L of the supernatant, transfer only the supernatant to another plate for measurement.

③ Preparation of dilution series of standard solution

1) Preparation of diluted buffer

Take 900 μ L each from the NAC Buffer (pH 8.0) and NAL Buffer (pH 10.2), and dispense each to another sample track. Add 60 μ L of water, 6 μ L of TFA and 234 μ L of acetonitrile to each. (hereinafter referred to as Diluted NAC Buffer, Diluted NAL Buffer)

2) Preparation of NAC and NAL standard solution

Take 300 μ L each from the prepared NAC solution and NAL solution, and dispense each to another sample track. Add 20 μ L of water, 2 μ L of TFA and 78 μ L of acetonitrile to each.

(hereinafter referred to as NAC standard solution, NAL standard solution)

3) Next, prepare the dilution series as follows.

	Std.1	Std.2	Std.3	Std.4	Std.5	Std.6	Std.7
Diluted NAC Buffer/ Diluted NAL Buffer	—	150 μ L	150 μ L	150 μ L	150 μ L	150 μ L	150 μ L
NAC/NAL Standard Solution	150 μ L	150 μ L	—	—	—	—	—
Dilution series operation ①	—	—	Std.2 at 150 μ L	—	—	—	—
Dilution series operation ②	—	—	—	Std.3 at 150 μ L	—	—	—
Dilution series operation ③	—	—	—	—	Std.4 at 150 μ L	—	—
Dilution series operation ④	—	—	—	—	—	Std.5 at 150 μ L	—
Final concentration	5.0 μ M	2.5 μ M	1.25 μ M	0.625 μ M	0.313 μ M	0.156 μ M	0

【Assay and data analysis】

After the reaction is completed, perform an analysis with HPLC.

Details of the conditions of analysis and method of analysis are available at the specialized webpage below.

ADRA kit specialized webpage

<https://labchem-wako.fujifilm.com/jp/download/adra/kit/index.html>

【Caution】

<Caution for measurement>

- Store the reagents in the specified conditions and do not use expired items.
- Use the kit as soon as possible after opening.
- Do not use the container and accessories in this kit for any other purposes.

<Caution for disposal>

- In case of disposal, follow the applicable effluent standards and dispose in a proper manner.

【Storage and expiration date】

- Storage: Keep in the refrigerator and away from the light (do not freeze).
- Expiration date: See the product label.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://www.wako-chem.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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