

Room-  
Temperature  
Storage

Refrigerated  
Storage

For short storage such as cell transportation!

# Cellstor-S, Cellstor-W

Cellstor-S is a cell suspension and preservation solution, and Cellstor-W is a cell wash and preservation solution manufactured by Otsuka Pharmaceutical Factory, Inc. They are used for refrigerated and room-temperature storage. They do not contain any components of human or animal origin. They are reagents for research use and are not intended for medical use in humans or animals.



Otsuka Pharmaceutical Factory, Inc.

## Features

- Store cells at room or refrigerated temperature
- Human and animal origin-free

## Usage

Extract the solution with a needle and a syringe. The product bags are not designed for cell storage. When storing cells, use a container suitable for cell storage.

### Differences in Usage of Cellstor-S and Cellstor-W

Since Cellstor-S is formulated with dextran 40, it is not suitable for washing cells when performing centrifugation. In such cases, using "Cellstor-W" solution, which does not contain dextran 40, is recommended.

## Cellstor-S    Cellstor-W

Contains 5% dextran 40  
(To suppress cell sedimentation)

Dextran 40 free



Unsuitable for centrifugation    Suitable for centrifugation



## Composition


Ingredients and their quantities in a bag (250 mL)

| Ingredients              | Cellstor-S         | Cellstor-W |
|--------------------------|--------------------|------------|
| Dextran 40               | 12.5 g (5%)        | -          |
| Trehalose Hydrate        | 8.29 g (3%)        |            |
| Calcium Chloride Hydrate | 0.05 g (0.02%)     |            |
| Potassium Chloride       | 0.075 g (0.03%)    |            |
| Sodium Chloride          | 1.5 g (0.6%)       |            |
| Sodium L-Lactate         | 0.775 g (0.3%)     |            |
| pH Adjuster              | Appropriate Amount |            |
| Water for Injection      | Appropriate Amount |            |

| Product Number | Product Name | Package Size | Storage Condition |
|----------------|--------------|--------------|-------------------|
| 637-46391      | Cellstor-S   | 250 mL       | R.T.              |
| 630-46401      | Cellstor-W   | 250 mL       | R.T.              |


## Example of Use 1 ~Temporal Changes in Cell Viability~

Changes in the cell viability of human adipose-derived mesenchymal stromal cells (hADSCs) in Cellstor-S, Cellstor-W or various solutions after storage at 5°C or 25°C



### Temporal changes in cell viability (25°C)

#### Room-Temperature Storage



|                     |                          |
|---------------------|--------------------------|
| Cell species        | hADSCs                   |
| Cell concentration  | $5 \times 10^5$ cells/mL |
| Storage temperature | 25°C                     |


#### Cell viability after 24-hour storage

|            |       |
|------------|-------|
| Cellstor-S | 98.1% |
| Cellstor-W | 95.6% |

Mean±SD (n=4), \*, p<0.05 Dunnett's test vs lactated Ringer's solution

### Temporal changes in cell viability (5°C)

#### Refrigerated Storage



|                     |                          |
|---------------------|--------------------------|
| Cell species        | hADSCs                   |
| Cell concentration  | $5 \times 10^5$ cells/mL |
| Storage temperature | 5°C                      |

#### Cell viability after 24-hour storage


|            |       |
|------------|-------|
| Cellstor-S | 95.1% |
| Cellstor-W | 93.5% |

Mean±SD (n=4),  
 \*\*\*, p<0.001 Dunnett's test, preservation solution A, B, C vs Cellstor-S  
 †; p<0.05, †††; p<0.001 Dunnett's test, preservation solution A, B, C vs Cellstor-W


**Maintained high viability up to at least 24 hours**

## Example of Use 2 ~Colony Formation Ratio~

Colony-forming capacity of human adipose tissue-derived mesenchymal stromal cells after 6-hour or 24-hour storage at 5°C or 25°C in Cellstor-S or Cellstor-W




Colony formation ratio before and after storage in Cellstor-S



Mean±SD (n=4)

There was no statistical significance ( $\alpha=0.05$ ) Dunnett's test vs Before preservation

Colony formation ratio before and after storage in Cellstor-W



Mean±SD (n=4)


There was no statistical significance ( $\alpha=0.05$ ) Dunnett's test vs Before preservation

**Maintained colony formation ratio up to at least 24 hours**


## Example of Use 3 ~Confirmation of Adipogenic Differentiation

Representative images of the adipogenesis differentiation assay with hADSCs ( $5 \times 10^5$  cells/mL) preserved in Cellstor-S and Cellstor-W after 24-hour storage at 5°C or 25°C. Adipogenic differentiation was induced and evaluated by Oil Red O staining.

Adipogenic differentiation before and after storage in Cellstor-S




Adipogenic differentiation before and after storage in Cellstor-W




## Example of Use 4 ~Confirmation of Osteogenesis Differentiation Ability~

Representative images of the osteogenesis differentiation assay with hADSCs ( $5 \times 10^5$  cells/mL) preserved in Cellstor-S and Cellstor-W after 24-hour storage at 5°C or 25°C. Osteogenic differentiation was induced and evaluated with an alkaline phosphatase staining kit and a calcified nodule staining kit.



## Example of Use 5 ~Confirmation of Cell Suspension Ability~

The percentage of cells in the supernatant after 1 hour of settling in normal saline, lactated Ringer's solution and Cellstor-S (The cells in the supernatant immediately after suspension were treated as 100%.)



**Cellstor-S prevents cell sedimentation for at least 1 hour**

This leaflet is based on information as of March 2021.

Listed products are intended for laboratory research use only, and not to be used for drug, food or human use. / Please visit FUJIFILM Wako Laboratory Chemicals site: <https://labchem-wako.fujifilm.com/> / This leaflet may contain products that cannot be exported to your country due to regulations. / Bulk quote requests for some products are welcomed. Please contact us.

**FUJIFILM Wako Laboratory Chemicals site**  
<https://labchem-wako.fujifilm.com>

**FUJIFILM Wako Pure Chemical Corporation**  
 1-2, Doshomachi 3-Chome, Chuo-ku, Osaka 540-8605, Japan  
 Tel: +81 6 6203 3741 Fax: +81 6 6203 1999  
 ffwk-cservise@fujifilm.com

**FUJIFILM Wako Chemicals U.S.A. Corporation**  
 1600 Bellwood Road, Richmond, VA 23237, U.S.A.  
 Toll-Free (U.S. only): +1 877 714 1920  
 Tel: +1 804 271 7677 Fax: +1 804 271 7791  
 wkuslabchem@fujifilm.com

**FUJIFILM Wako Chemicals Europe GmbH**  
 Fuggerstr 12, 41468 Neuss, Germany  
 Tel: +49 2131 311 0 Fax: +49 2131 311 100  
 labchem\_wkeu@fujifilm.com

**FUJIFILM Wako Chemicals (Hong Kong) Limited**  
 Room 1111, 11/F, International Trade Centre, 11-19 Sha Tsui Road,  
 Tsuen Wan, N.T., Hong Kong  
 Tel: +852-2799-9019 Fax: +852-2799-9808  
 wkhk.info@fujifilm.com

**FUJIFILM Wako (Guangzhou) Trading Corporation**  
 Room 3003, 30/F., Dong Shan Plaza 69, Xian Lie Zhong Road,  
 Guangzhou, 510095, China  
 Tel: +86-20-8732-6381(Guangzhou) Tel: +86-21-6288-4751(Shanghai)  
 Tel: +86-10-6413-6388(Beijing)  
 wkgz.info@fujifilm.com