Development of a novel highly sensitive oxytocin ELISA


*Shiro Sugino, Masaaki Kojima, Fumie Akutsu, Takahiro Nishibu FUJIFILM Wako Pure Chemical Corporation



292-84401 Oxytocin ELISA kit *Wako*

Background

Oxytocin is a peptide hormone comprising nine amino acids. It is produced in the hypothalamus and released mainly from the posterior pituitary gland following childbirth and during lactation, and promotes uterine contraction and milk production. It also has stress-relieving, anxiolytic, and fear-reducing effects and plays a role in the development of maternal behavior. Due to these properties, oxytocin is often referred to as the "happy hormone" or "love hormone." Recently, oxytocin has attracted attention in treating mental disorders such as depression and autism and in the development of functional food materials.



Lactation
Uterine contraction
Milk production

Stress-relieving
Anxiolytic
Fear-reducing
Social behavior
Anti-obesity


Commercially available ELISA kits require a large volume of sample and complicated pretreatment using a C18 column and organic solvent.




We developed a novel ELISA for oxytocin needing only simple pretreatment and a small sample.

Assay principal

We constructed an oxytocin competition ELISA with high sensitivity (4 pg/mL) using biotin-labeled, streptavidin-conjugated HRP and luminescent substrate.





Sample pretreatment

Oxytocin binds to proteins and other molecules in the sample, and these can affect measurement by inhibiting oxytocin-antibody interaction (i.e., interactions between oxytocin and antibodies). We investigated various methods for removing the impurities in the sample and found the following novel pretreatment method without extraction process using organic solvent.


Fujifilm Wako Sample (saliva, serum, plasma, urine) Add sample buffer 1 (acid) Mix, and incubate 5 min at R.T., and mix incubate 5 min at R.T. Add sample buffer 2 (gel) Mix, and incubate 5 min at R.T., and mix incubate 5 min at R.T., Centrifuge Use the supernatant as a sample for measurement.

Conventional method	
Sample (saliva, serum, plasma, urine)	
Add TFA, and mix Centrifuge	
Add the supernatant to C18 Column	
equilibrated with TFA and acetonitrile	
Wash with TFA, and elute with acetonitrile	
Depressurize and dry	
Dissolve and use the solution as a sample	

	Fujifilm Wako Conventional method		
Preparation time	30 min	2 h-overnight	
Sample volume	50-200 μL	250 μL-3 mL	
Recovery rate	90-120%	42-110%	


Spiked-recovery test


Our ELISA demonstrated good recovery rate (near 100%) using human saliva, serum, and urine-spiked oxytocin. On the other hand, the conventional method shows too low recovery rate in case of serum, and non-pretreated method shows too high recovery rate in the all samples.



Dilution linearity


Our ELISA demonstrated good dilution linearity using human plasma, serum, urine, and saliva, and mouse/rat plasma and serum. The representative data are shown.





Wako


Example of measurement -healthy human / mouse / rat samples-



Sample concentration method

■Method for 3-fold concentration

Sample (after preparation): 360 μ L



Reconstitute the dry samples with 120 μ L buffer

Measured by ELISA

Sample	ole No.	Before	After	After/before
Sample		Measurement (pg/mL)	Measurement (pg/mL)	After/before
Luman	1	28.4	80.5	2.83
Human Saliva	2	15.1	40.3	2.67
Jativa	3	20.5	55.2	2.69
Lluman	1	48.8	122	2.50
Human Serum	2	31.9	99.7	3.13
JCIUIII	3	33.8	95.6	2.83
	1	36.8	105	2.85
Human Urine	2	12.0	31.9	2.66
Office	3	44.5	151	3.39

3-fold concentration method induces about triple the measured value.

This method is useful for measuring samples with lower oxytocin concentration.

Conclusion

- We have developed a novel oxytocin ELISA method that can measure a variety of samples using a simple pretreatment method.
- Our method is expected to be a useful tool for various oxytocin studies.
- We are developing ELISA using monoclonal antibody aiming stable supply as a next step.

Disclosure of Conflict of Interest

Matters requiring disclosures of COI with regard to the poster are as follows: Employee of FUJIFILM Wako Pure Chemical Corporation