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1. Abstract


4. Evaluation of Fg in the novel culture method

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Outline of this study


- ◆ Improve gastrointestinal drug metabolism
- ◆ Improve *in vitro* Fg prediction accuracy
- ◆ Make F-hiSIEC™ a more useful tool for DMPK and TOX


1) Imakura, et al. *Biochem Biophys Res Commun.* 692:149356, 2024
 2) Shirai, et al. *Drug Metab Pharmacokinet.* 55:100994, 2024

2. Method

Overview of differentiation protocol


User culture protocol (trans well culture)

Difference between novel and conventional methods

	Cell	Medium		Culture method
		Seeding	Culture	
Novel ²⁾	F-hiSIEC™ Cell	F-hiSIEC™ Seeding	New	Day0-4: LL Day4-11: ALI
Conv.	Cell	Seeding	F-hiSIEC™ Culture	Day 0-11: LL


2) Shirai et al., *Drug Metab Pharmacokinet* 55: 100994 (2024)

Formula for calculating *in vitro* Fg³⁾

$$Fg = \frac{A \text{ to } B \text{ transport}}{(A \text{ to } B \text{ transport}) + (\text{Metabolite production})}$$


3) Michiba et al., *Drug Metab Dispos* 50: 204-213 (2022)


- Estimated Fg value of each compound were calculated by the equation described above.
- A to B transport represents the transcellular transport clearance in the apical-to-basolateral direction.
- Metabolic production was calculated by quantifying metabolite in apical, cell and basal compartment.

3. Gene expression, Barrier function

Gene expression

Barrier function


✓ The novel method increased the expression of small intestinal markers, transporters, and metabolizing enzymes.

✓ Barrier function was maintained in the novel method as in the conventional method.

5. Fg prediction of CYP3A substrates

Estimated Fg values for typical CYP3A4 substrates

Substrate	Method	A to B transport (μL/hr/well)		Metabolite production (μL/hr/well)		Estimated Fg value		<i>In vivo</i> human Fg value ⁴⁾⁵⁾
		Mean	SD	Mean	SD	Mean	SD	
Felodipine	Conventional	16.5	1.4	2.13	0.36	0.88	0.02	0.35
	Novel	7.22	1.3	10.1	1.70	0.42	0.04	
Midazolam	Conventional	60.6	8.7	7.43	0.54	0.89	0.02	0.55
	Novel	39.9	3.0	18.6	2.17	0.68	0.02	
Verapamil	Conventional	26.4	7.9	0.64	0.24	0.98	0.00	0.73
	Novel	30.0	1.6	4.71	0.62	0.86	0.02	
Sildenafil	Conventional	16.9	7.1	0.390	0.01	0.98	0.01	0.78
	Novel	17.7	0.4	1.97	0.20	0.90	0.01	

4) Verma et al., *J Med Chem* 53: 1098-1108 (2010) 5) Gertz et al., *Curr Drug Metab* 9: 785-795 (2011)

✓ The Fg values predicted by the novel method using F-hiSIEC™ showed good correlation with measured human values.

6. Potential for expansion into toxicity assessment

Various applications for gastrointestinal research have been developed in F-hiSIEC™.

- Intestinal toxicity evaluation** *Report at 48th and 49th JSOT
 - Evaluation of gastrointestinal cytotoxicity (e.g., 5-FU and its prodrugs)
 - Evaluation of differentiation abnormalities (e.g., increase in goblet cells with γ-secretase inhibitor)
 - Decreased barrier function in inflammatory conditions (e.g., addition of inflammatory cytokines)
- DMPK evaluation**
 - Fg prediction (Fg prediction for CYP3A substrate) ⇒ Novel culture method reported here
 - Fa prediction (P_{app} evaluation on various drugs)
 - Enzyme induction (CYP3A4 gene induction by RIF and VD3).
- Immune and inflammatory evaluation**
 - Evaluation of substance uptake via M cells
 - Evaluation of inflammation suppression by short chain fatty acids
- Co-culture with intestinal bacteria**
 - Co-culture with *B. fragilis*, *Lactobacillus*⁶⁾, and *Salmonella*⁷⁾ and evaluation of interactions
- Norovirus culture, disinfectant evaluation**

6) Sen, A. et al., *Front Microbiol.* 2023 Apr 13:14:1155438
 7) Fuka, Y. et al., *FEMS Microbiology Letters*, 2024, 371, fnae006

- ✓ F-hiSIEC™ could be used to assess toxicity through drug metabolism and to predict systemic exposure for toxicity assessment.
- ✓ Evaluation of drug toxicity via intestinal bacteria using F-hiSIEC™ is also expected.

Conclusion : The novel culture method improved the metabolic function of F-hiSIEC™ and the accuracy of Fg prediction of CYP3A substrate drugs. This evaluation system may be useful as a model for evaluating gastrointestinal absorption. In addition, the ability to accurately predict gastrointestinal metabolism and Fg may be applicable to toxicity assessment.